

MULTIGATE MamaCool™



Perineal cooling pads that provide safe and effective pain relief for mums

RE-ORDER
10-010
Qty.
10 pieces/box
12 boxes/carton

FEATURES

Cooling therapy

The MamaCool™ pad provides a safe, effective and convenient method of delivering cooling therapy to women experiencing perineal trauma post vaginal childbirth.

It's unique design provides a hygienic and comfortable pain management solution that supports mothers in their post-partum recovery.



Reduces perineal pain and swelling without medication

INNOVATION

Designed by a midwife

Throughout her 25 years as a midwife, Natalie was frustrated by the lack of a dignified, fit for purpose product to deliver relieving cold therapy to new mums in her care.

This led her to develop the MamaCool™ pad so that clinicians could easily follow cooling therapy protocols and women's recovery experiences could be improved.

- ✓ Provides therapeutic levels of pain relief without medication¹
- ✓ Hygienic, mess free and easy to use solution
- ✓ Doesn't interfere with breastfeeding



NATALIE W.
PRODUCT INVENTOR
& MIDWIFE

"Women need a comfortable option that reduces pain and swelling. Midwives want a product that is easy to use, hygienic and works."

¹ In a trial, MamaCool™ was shown to reduce perineal temperatures by an average of 11.7 degrees Celsius, thus achieving therapeutic levels of cooling.


For more information please contact your Business Development Manager:
(02) 8718 2888 or visit multigate.com.au

DESIGN

Maximum comfort

MamaCool™ was created with the woman's comfort in mind. It's distinctive holes allow post-birth blood to flow through to the maternity pad underneath to keep the vagina and

perineum dry. The soft and smooth semi-absorbent outer layer is gentle on stitches and the shape of the pad contours to the woman's body.



SUPPORT


Instructions for use

- 1** ❄️ Store MamaCool™ pads in the freezer ready for use (minimum 1 hour). The pad should be firm.
- 2** ⌚ Remove a MamaCool™ pad from the freezer and allow it to soften slightly. This will take approximately 5-10 minutes.
- 3** 📎 Place MamaCool™ on top of maternity pad* and wear together for 10-20 minutes (or longer if desired). The orientation and position of the MamaCool™ pad can be adjusted to align with the area where cooling relief is required.

- 4** 🗑️ Remove the MamaCool™ pad and dispose.
 - 5** 🔄 Reapply a new MamaCool™ pad every 2-4 hours for the first 72 hours after birth and as often as desired thereafter.
- !** Store in a dry, cool and well-ventilated place if not storing in the freezer. Avoid direct sunlight.
- ⌚** MamaCool™ pads are a single use product. Dispose of pads once used.

*Maternity pad not included.

Tips



Scan the QR code for more information about caring for your perineum after child birth or visit www.multigate.com.au/qr06

Warning: Cold therapy should not be used by people with known or suspected circulatory problems unless under supervision of a doctor. If the perineal cooling pad is damaged and the contents contact skin, rinse off with water. Seek medical advice if irritation occurs and persists.

IP Protected Product: Intellectual property rights such as patents, patent applications, registered designs or trademarks apply to this product.