


Clip Tourniquet

Tighten infection control.

Our slow & quick release clip tourniquet provides an economic single use or single patient use solution to venous access to suit your unique hospital infection control protocols.



DESIGN

With growing evidence that tourniquets are a source of MRO transmission and cross-infection, our disposable clip tourniquet has been designed to complement meticulous sterile technique.

- ✓ Easy to use
- ✓ Slow and quick release
- ✓ Reduces cross-contamination
- ✓ Soft elastic for patient comfort
- ✓ Latex free



RANGE

Made from soft elastic, our clip tourniquet is comfortable when applied and won't damage fragile skin. Available as either a single use or single patient use item in adult and paediatric size.



The **single use** model is identifiable by its white clip which contains the “do not reuse” symbol.



The **single patient use** model has a teal (adult) or blue (paediatric) clip. The white band enables easier detection of soiling or contamination. Should this occur, please discard and replace with a new tourniquet.


CLIP TOURNIQUETS - Non-sterile

REORDER	DESCRIPTION	PIECE PER BOX	BOXES PER CARTON	TOTAL PIECES
DP 3509	Clip Tourniquet - Adult (Single Use)	50pcs	10boxes	500pcs
25-116NS	Clip Tourniquet - Adult (Single Patient Use)	50pcs	10boxes	500pcs
DP 3513	Clip Tourniquet - Paediatric (Single Use)	50pcs	10boxes	500pcs
25-117NS	Clip Tourniquet - Paediatric (Single Patient Use)	50pcs	10boxes	500pcs


Clip Tourniquet

Instructions for use


HANDLING




Wrap the tourniquet around the limb 7.5cm to 10cm above puncture site and fasten by inserting the clasp into the top slot of the buckle.



To tighten, pull the tourniquet band away from the buckle until the desired pressure is achieved.



For a slow release or to loosen the tourniquet, press the bottom ridged button and slightly lift the buckle.



For a quick release, press the top smooth button on the buckle.

CAUTIONS / PRECAUTIONS


Do not apply high pressures that may impair arterial blood flow to the extremity.

Release tourniquet as soon as blood starts to flow into the first tube.

Too long and too high stasis may contribute to haemolysis.



Avoid exposure to direct sunlight.



Store in a cool, dry and clean place.



Latex Free.