

Clip Tourniquet

Tighten infection control.

With growing evidence that tourniquets are a source of MRO transmission and cross-infection, our single use or single patient use Clip Tourniquets provide a superior solution to venous access that complements meticulous sterile technique.



DESIGN

Our slow & quick release Clip Tourniquet provides a familiar design that is easy to use and control. It allows for one-handed slow or quick release and strap detachment. Made from soft elastic, our clip tourniquet is comfortable when applied and won't damage fragile skin. Available in adult and paediatric size.

- ✓ Easy to use
- ✓ Slow and quick release
- ✓ Reduces cross-contamination
- ✓ Soft elastic for patient comfort
- ✓ Latex free



RANGE

Available as either a single use or single patient use item in adult and paediatric size.



The **single use** model contains the “do not reuse” symbol on the clip and a teal or light blue band.



The **single patient use** model contains the words “Single Patient Use” on the clip and a dark blue band.


CLIP TOURNIQUETS - Non-sterile

REORDER	DESCRIPTION	PIECE PER BOX	BOXES PER CARTON	TOTAL PIECES
DP 3509	Clip Tourniquet - Adult (Single Use)	50pcs	10boxes	500pcs
25-116NS	Clip Tourniquet - Adult (Single Patient Use)	50pcs	10boxes	500pcs
DP 3513	Clip Tourniquet - Paediatric (Single Use)	50pcs	10boxes	500pcs
25-117NS	Clip Tourniquet - Paediatric (Single Patient Use)	50pcs	10boxes	500pcs


Clip Tourniquet

Instructions for use


HANDLING




Wrap the tourniquet around the limb 7.5cm to 10cm above puncture site and fasten by inserting the clasp into the top slot of the buckle.



To tighten, pull the tourniquet band away from the buckle until the desired pressure is achieved.



For a slow release or to loosen the tourniquet, press the bottom ridged button and slightly lift the buckle.



For a quick release, press the top smooth button on the buckle.

CAUTIONS / PRECAUTIONS


Do not apply high pressures that may impair arterial blood flow to the extremity.

Release tourniquet as soon as blood starts to flow into the first tube.

Too long and too high stasis may contribute to haemolysis.



Avoid exposure to direct sunlight.



Store in a cool, dry and clean place.



Latex Free.