



Abdominal sponges & swabs

DESIGNED TO STAND OUT

Ensuring patient safety is paramount to us. That's why our Fluro-Tec® range is designed to stand out and complement your processes to help lower the risk of retained sponges and swabs.

FEATURES AND BENEFITS

- ✓ Unique, patented fluorescent thread
- ✓ Contrast against blood improves intra-operative detection
- ✓ Thread glows when exposed to UV light
- ✓ Highly absorbent
- ✓ Low linting and no loose threads



Intra-operative detection — Fluro-Tec® sponges and swabs contain a fluorescent, non-absorbent green-yellow thread to make a blood soaked Fluro-Tec® 9.5 times more visible¹ to the naked eye and even more detectable under UV light.




BETTER FOR THE PLANET

Our composite sponges can reduce greenhouse gas emissions by up to 25%² and they feel and absorb the same as our regular cotton sponge.



IMPROVED
PROTECTION

Extra protection against 'never events' — our unique fluorescent thread ensures our sponges and swabs are visible under UV light and detectable by the naked eye even when soaked in blood.



X-RAY DETECTABLE 6 PLY ABDOMINAL SPONGE

RE-ORDER	PIECES (PCS)	TYPE	SIZE (CM)	STERILITY	UNITS PER INNER (PCS)	UNITS PER CARTON
88-987	5	Composite	22.5 x 22.5	Sterile	-	60 Pks
88-988	5	Composite	30 x 30	Sterile	-	50 Pks
88-989	5	Composite	50 x 30	Sterile	-	30 Pks
88-987NS	1	Composite	22.5 x 22.5	Non-Sterile	100	8 Bags
88-988NS	1	Composite	30 x 30	Non-Sterile	100	4 Bags
88-989NS	1	Composite	30 x 50	Non-Sterile	100	3 Bags

X-RAY DETECTABLE 8 PLY SWABS

RE-ORDER	PIECES (PCS)	TYPE	SIZE (CM)	STERILITY	UNITS PER INNER (PCS)	UNITS PER CARTON
88-683	5	Standard	10 x 10	Sterile	-	250 Pks
88-684	5	Premium	10 x 10	Sterile	-	250 Pks
88-684NS	5	Premium	10 x 10	Non-Sterile	-	800 Pks

A BRIGHT IDEA

Developed in Australia in close consultation with clinicians, this product reflects real-world needs translated into smarter design.

Disclaimers: 1. Luminance contrast performance of Fluro-Tec® compared to a standard swab or sponge, independently verified under a control study "Multigate Fluro-Tec® Luminance Contrast Assessment", Ergonomie Pty Ltd, 04/07/2019. 2. Up to 25% reduction in Greenhouse Gas (GHG) emissions compared to an equivalent Multigate cotton sponge, based on End of Life (EOL) clinical waste disposal method 100% disinfect and landfill being used. Savings from using EOL 100% incineration is 17%.

FOR MORE INFORMATION

multigate.com.au or call 1800 023 420 (Australia)
multigate.co.nz or call 0800 880 013 (New Zealand)