

# MULTIGATE SafeCount® Needle Counters



## NEEDLE COUNTERS → SAFECOUNT® TRANSITION TABLE (STERILE)

Needle Counters

PRODUCT DESCRIPTION	Needle Counter 20CT Foam Strip Clear Lid	Needle Counter 30CT Foam Strip Clear Lid	Needle Counter 20CT Foam Strip Clear Lid	Needle Counter Magnetic	Needle Counter 30CT Magnetic	Needle Counter 15CT Foam Strip Magnetic	Needle Counter 42CT Foam Strip	Needle Counter 20CT Foam Magnetic	Needle Counter 30CT Foam Magnetic	Needle Counter 40CT Foam Magnetic
SIZE WHEN CLOSED	140 x 72 x 15mm	140 x 72 x 15mm	145 x 52 x 34mm	141 x 72 x 14mm	127 x 52 x 16mm	127 x 53 x 16mm	140 x 72 x 15mm	145 x 52 x 34mm	140 x 72 x 15mm	153 x 91 x 37mm
PIECE WEIGHT	60g	58g	61g	85g	61g	53g	65g	76g	77g	122g
CODE	19-233	19-237	19-241	19-231	19-236	19-232	19-240	19-234	19-238	19-239
GHG D&L*	351 gC02e	346 gC02e	353 gC02e	453 gC02e	332 gC02e	300 gC02e	373 gC02e	426 gC02e	433 gC02e	712 gC02e
GHG I*	487 gC02e	479 gC02e	492 gC02e	650 gC02e	472 gC02e	422 gC02e	521 gC02e	600 gC02e	609 gC02e	992 gC02e



## TRANSITIONING TO

SafeCount®

PRODUCT DESCRIPTION	Needle Counter 20 Count 20 Capacity Hybrid with Clear Lid	Needle counter 20 Count 40 Capacity Double Magnetic	Needle Counter 20 Count 40 Capacity Magnetic and Hybrid	Needle Counter 20 Count 40 Capacity Magnetic and Foam	Needle Counter 45 Count 90 Capacity Magnetic and Foam
SIZE WHEN CLOSED	127 x 59 x 16mm	127 x 59 x 16mm	127 x 59 x 16mm	127 x 59 x 36mm	153 x 73 x 36mm
PIECE WEIGHT	38g	64g	56g	70g	100g
CODE	CEA-0014S	CEA-0015S	CEA-0016S	CEA-0010S	CEA-0011S
GHG D&L*	270 gC02e	388 gC02e	337 gC02e	424 gC02e	684 gC02e
GHG I*	369 gC02e	534 gC02e	464 gC02e	584 gC02e	914 gC02e
GHG CHANGE UP TO:	-24.32% [20CT Foam Strip Clear Lid 19-233] -22.94% [30CT Foam Strip Clear Lid 19-237] -25.08% [20CT Foam Strip Clear Lid 19-241]	-17.81% [Magnetic 19-231] 16.88% [30CT Magnetic 19-236]	-12.32% [15CT Foam Strip Magnetic 19-232] -11.05% [42CT Foam Strip 19-240]	-2.61% [20CT Foam Magnetic 19-234]	58.15% [30CT Foam Magnetic 19-238] -7.87% [40CT Foam Magnetic 19-239]

TABLE KEY: (GHG) Greenhouse Gas Emissions, (D&L) Disinfect & Landfill, (I) Incineration.

### Benefits of needle counter transition

SafeCount® needle counters provide better management of sutures, blades, sharps and handles.




#### Sharps safety

Slots in foam allow for easy insertion and removal of sutures.



#### Easy blade lift

Undulating magnetic strip assists with blade retrieval from magnet.



#### Instrument organisation

Notches at one end of the unit can be used to rest instruments.



As part of our commitment to better product stewardship, all SafeCount® needle counters are supplied climate neutral<sup>1</sup>.



We use 100% recycled paper when we print this item and we offset its greenhouse gas emissions to make it climate neutral

For more information please contact your local Business Development Manager on (02) 8718 2888 or visit [multigate.com.au](https://multigate.com.au)

\*Estimated total GHG emissions figures represent the cradle-to-grave life cycle impact of the product for the nominated End Of Life disposal method. As all figures were created using the same LCA model, methodology, data sources and assumptions, an accurate apple-for-apple comparison of products is possible. LCA figures from different sources should not be compared as modelling differences could result in an unreliable comparison. 1. Multigate purchases carbon credits via Carbon Neutral in approved, registered projects which have been independently assessed and verified as creating additional resources that help combat climate change using avoidance or sequestration means. LCA data contained is protected by copyright and cannot be reproduced or adapted with prior written permission from Multigate. IP Protected Product see website for details.