



# Surgical Face Mask (Level 2)



## DESCRIPTION

Surgical masks, also known as medical procedure masks, are used to prevent large droplets of bodily fluids that may contain viruses or bacteria from escaping the wearer's nose and mouth. It helps to reduce the chance of spreading infectious airborne

viruses and bacteria to others. It may also help block large droplets, splashes, sprays or splatter that may contain viruses and bacteria from others, keeping it from reaching the wearer's nose and mouth.

## FEATURES AND BENEFITS

- ✓ AS4381:2015 compliant
- ✓ Bacterial Filtration Efficiency (BFE)  $\geq$  98%
- ✓ Fluid resistant - 120mmHg (Level 2)
- ✓ 3 ply Construction
- ✓ Single use

## PRODUCT DATA

- 17.5cm x 9.5cm
- Non-woven melt-blown material
- Non-latex elastic band
- Expansion pleats (covers from nose to chin)
- Bendable nose clip and ear loops for securement

## FITTING INSTRUCTIONS

1



Place surgical mask over your face with the metal strip sitting on top of your nose and position both elastic bands behind your ears.

2



Mould the metal strip to the shape of your nose to ensure it fits your face comfortably.

3



Adjust the mask to ensure it covers your nose and chin.

## Personal Protection Equipments

RE-ORDER	DESCRIPTION	SIZE (CM)	STERILE / NON-STERILE	PIECES PER INNER BOX	INNER PER CARTON	TOTAL PIECES PER CARTON
54-420	Surgical Face Mask with ear loops	17.5cm x 9.5cm	Non-Sterile	50pcs	40 Boxes	2000pcs