

All Purpose Towels

An economic alternative

A versatile towel that is highly absorbent with fast fluid capture rate. Can be used wet or dry.

Features and benefits

High absorbency with quick saturation rate makes our All Purpose Towels suitable for a wide range of cleaning applications.

- ✓ Absorbs up to 8.5 times its weight in fluid
- ✓ Fast fluid capture
- ✓ Superior wet and dry strength
- ✓ Low linting
- ✓ Latex free

MediClean™ - Product Specifications (4plys folded)

MATERIAL/ CONSTRUCTION	WEIGHT OF TOWEL	WATER ABSORBENCY	SINKING TIME	TENSILE STRENGTH: Machine Direction	TENSILE STRENGTH: Cross Direction
50% Viscose 50% Polyester	45gsm	8.5 g/g	< 1.3 sec	98 N	28 N



MediClean™ - All Purpose Towel

RE-ORDER	DESCRIPTION	SIZE (CM)	STERILE / NON-STERILE	PIECES PER DISPENSER	PIECES PER CARTON
57-885	All Purpose Towels Small	35 x 30	Non-Sterile	100pcs	1200pcs
57-886	All Purpose Towels Large	35 x 60	Non-Sterile	100pcs	1200pcs