

Catheter Pack

06-687

Enhancing performance at every step.

Our Catheter Packs and in-service support ensure an efficient, managed process designed to support and enhance catheterisation.

PERFORMANCE

Improve Sterile Protocols.

Managing sterile protocols is paramount. All our Catheter Packs are supported by an in-service education and training support program to ensure sterile protocols, patient comfort and product performance are maximised.

Reduce Waste.

Our Catheter Packs support standardised and compliant practice and optimise the use of components for your catheterisation procedure.

Pack Contents Catheter Pack (Sterile)

RE-ORDER: 06-687 CARTON QUANTITY: 25 PACKS

QTY	DESCRIPTION	REF.
1	Sterile Field (90cm x 90cm)	A
2	Kidney dish 700mL – natural	H
1	Graduated Galipot 150mL	G
5	XRD Gauze Swabs (10cm x 10cm)	I
1	Towel (35cm x 60cm)	C
1	Drape (95cm x 50cm) with fenestration (20cm x 8cm)	B
2	Forceps (plastic) 13cm – Blue	D
1	Lubricating jelly 2.7gm	F
1	Injectable Water 10mL H ₂ O Steritube	E

- ✓ Optimal component usage
- ✓ Supporting sterile protocols
- ✓ Improved workflow
- ✓ Improved patient safety

