

# MULTIGATE

## Claripose®



### Adhesive Island Dressings

#### Island Dressings Range

Our Claripose® island dressings are ideal for use as a post-operative wound dressing. The “all-in-one” dressing consists of a non-adherent absorbent pad attached to a transparent, waterproof adhesive backing which is breathable yet impermeable to liquids, bacteria and viruses. The flexible backing helps the dressing conform to the body, making it suitable for use on high movement body areas.



- ✓ Sterile, flexible & breathable
- ✓ Waterproof and provides bacterial & viral barrier
- ✓ Conformable around difficult anatomical contours
- ✓ Secure and reliable adhesive, gentle to the skin
- ✓ Easy application and removal
- ✓ Latex free

#### Claripose® Island Dressings (Peel pack of 1)

RE-ORDER	DESCRIPTION	SIZE	STERILE / NON-STERILE	PACK PER BOX	BOX PER CARTON	PACK PER CARTON
30-651	Claripose® Adhesive Island Dressing	6.5cm x 5cm	Sterile	100	6	600pks
30-652	Claripose® Adhesive Island Dressing	9.5cm x 8.5cm	Sterile	20	10	200pks
30-653	Claripose® Adhesive Island Dressing	12cm x 10cm	Sterile	10	10	100pks
30-654	Claripose® Adhesive Island Dressing	15.5cm x 8.5cm	Sterile	20	10	200pks
30-658	Claripose® Adhesive Island Dressing	34cm x 10cm	Sterile	20	8	160pks

For more information please contact your local Business Development Manager:

(02) 8718 2888 or visit [multigate.com.au](http://multigate.com.au)