

MULTIGATE

Compose®



Adhesive Island Dressings

Island Dressings Range

Our Compose® Island Dressing is a soft, comfortable, self-adhesive dressing designed for use with all wound types. The dressing conforms easily to anatomical contours and is perfect for post-operative wounds. It consists of an Interpose® non-adherent pad and a fluid repellent, non-woven backing layer that allows the movement of air. It also features a low irritant adhesive border that provides dressing security and allows for easy removal.



- ✓ Sterile, highly absorbent, non-adherent pad
- ✓ Soft, flexible cover that is breathable and conforms to anatomical contours
- ✓ Secure and reliable adhesive, gentle on the skin

- ✓ Easy application and removal
- ✓ Patient comfort
- ✓ Latex free

Compose® Island Dressings (Peel pack of 1)

RE-ORDER	DESCRIPTION	SIZE	STERILE / NON-STERILE	PACK PER BOX	BOX PER CARTON	PACK PER CARTON
30-611	Compose® Adhesive Island Dressing	7.5cm x 5cm	Sterile	50	8	400pks
30-612	Compose® Adhesive Island Dressing	8.25cm x 6cm	Sterile	50	8	400pks
30-613	Compose® Adhesive Island Dressing	12cm x 8.25cm	Sterile	20	8	160pks
30-614	Compose® Adhesive Island Dressing	18cm x 8.25cm	Sterile	20	8	160pks
30-615	Compose® Adhesive Island Dressing	25cm x 8.25cm	Sterile	20	8	160pks

For more information please contact your local Business Development Manager:
(02) 8718 2888 or visit multigate.com.au