

MULTIGATE SmartTab[®] Dressing Packs



Keeping tabs on infection.

Powered by the SmartTab[®] Aseptic Field, Multigate Dressing Packs have been assessed by the Association for Safe Aseptic Practice as the first in Australia to officially support ANTT[®] (Aseptic Non Touch Technique).

INNOVATION

Supporting aseptic practice.

SmartTab[®] Dressing Packs help to protect Key-Parts and Key-Sites in any procedure.

- ✓ Pack folding supports a safe and efficient opening
- ✓ Tabs provide a simple and effective way to open and position the Critical Aseptic Field
- ✓ Two-colour system clearly denotes aseptic surfaces

MULTIGATE
SmartTab[®]
ANTT
COMPLIANT
DRESSING PACKS




DESIGN

Improved process.

SmartTab[®] Dressing Packs support an easy four-step preparation of the Aseptic Field that helps caregivers to consistently follow ANTT[®] principles.



Visit: multigate.com.au/smarttab



FEATURES

Confidence.

Our forceps are designed to prevent lateral movement of the arms, preventing anything clamped from slipping free.

Security.

Packs are folded and secured with a tear-open sticker, reducing the risk of unintentional unravelling of the Pack.

Resources.

Formal ANTT[®] compliance better connects our customers to best practice resources prepared by The Association for Safe Aseptic Practice.

For more information please contact your Business Development Manager:
(02) 8718 2888 or visit multigate.com.au

**SMALL CHANGE
= big IMPACT**

MULTIGATE

SmartTab[®]

Dressing Packs



RE-ORDER

00-001

Supporting aseptic practice

- ✓ Pack folding supports a safe and efficient opening
- ✓ Tabs provide a simple and effective way to open and position the Critical Aseptic Field
- ✓ Two-colour system clearly denotes aseptic surfaces

SCAN THE QR CODE BELOW TO WATCH AN INSTRUCTION VIDEO



MULTIGATE
SmartTab[®]
ANTT[®]
COMPLIANT
DRESSING PACKS



Dressing Packs Tear Pack 160 Packs/Carton

DESCRIPTION	* KEY-PARTS
1 x Aseptic field * 1 x Dressing tray * 1 x Paper towel 2 x Non-woven swabs * 5 x Non-woven Balls * 1 x Yellow forcep * 2 x Blue forceps *	Key-Parts are those parts of the equipment that come into contact with any open wounds, insertion and puncture sites, liquid infusion or with other active Key-Parts connected to the patient via a medical device.




Flared sides are designed to provide stability.



Sloping grid aids the removal of excess fluid.



High walls minimise chance of spills and accommodate generous fluid volumes.



For more information please contact your Business Development Manager:

(02) 8718 2888 or visit multigate.com.au

**SMALL CHANGE
= big IMPACT**