

MULTIGATE SmartTab[®] Dressing Packs



Keeping tabs on infection.

Powered by the SmartTab[®] Aseptic Field, Multigate Dressing Packs have been assessed by the Association for Safe Aseptic Practice as the first in Australia to officially support ANTT[®] (Aseptic Non Touch Technique).

INNOVATION

Supporting aseptic practice.

SmartTab[®] Dressing Packs help to protect Key-Parts and Key-Sites in any procedure.

- ✓ Pack folding supports a safe and efficient opening
- ✓ Tabs provide a simple and effective way to open and position the Critical Aseptic Field
- ✓ Two-colour system clearly denotes aseptic surfaces

MULTIGATE
SmartTab[®]
ANTT
COMPLIANT
DRESSING PACKS




DESIGN

Improved process.

SmartTab[®] Dressing Packs support an easy four-step preparation of the Aseptic Field that helps caregivers to consistently follow ANTT[®] principles.



Visit: multigate.com.au/smarttab



FEATURES

Confidence.

Our forceps are designed to prevent lateral movement of the arms, preventing anything clamped from slipping free.

Security.

Packs are folded and secured with a tear-open sticker, reducing the risk of unintentional unravelling of the Pack.

Resources.

Formal ANTT[®] compliance better connects our customers to best practice resources prepared by The Association for Safe Aseptic Practice.

For more information please contact your Business Development Manager:

(02) 8718 2888 or visit multigate.com.au

DESIGN

Supporting your needs.

Our Dressing Pack configurations support a range of clinical preferences.



Managing fluids.

Our two dressing trays are designed to avoid spills and improve the management of fluids.



Flared sides are designed to provide stability.



Sloping grid aids the removal of excess fluid.



High walls minimise chance of spills and accommodate generous fluid volumes.



Tray 1



Tray 2

Tray 1: Dressing Packs Tear Packs		
RE-ORDER	DESCRIPTION	PER CARTON
08-559NP	1 aseptic field, 1 dressing tray, 1 paper towel, 5 non-woven swabs, 1 yellow forcep, 2 blue forceps	160
00-001	1 aseptic field, 1 dressing tray, 1 paper towel, 2 non-woven swabs, 5 non-woven balls, 1 yellow forcep, 2 blue forcep	160
08-888NP	1 aseptic field, 1 dressing tray, 1 paper towel, 2 gauze swabs, 5 cotton balls, 1 yellow forcep, 2 blue forcep	160
08-331P	1 aseptic field, 1 dressing tray, 1 paper towel, 6 non-woven swabs, 1 yellow forcep, 2 blue forceps	160
00-008	1 aseptic field, 1 dressing tray, 1 paper towel, 5 non-woven swabs, 1 yellow forcep, 2 blue forceps	160
07-888P	1 aseptic field, 1 dressing tray, 1 paper towel, 3 gauze swabs, 6 cotton balls, 1 yellow forcep, 2 blue forceps	120
00-001D	1 aseptic field, 1 dressing tray, 1 paper towel, 2 non-woven swabs, 5 non-woven balls, 1 yellow forcep, 2 blue forceps	90

Tray 1: Dressing Packs Peel Packs		
RE-ORDER	DESCRIPTION	PER CARTON
00-000P	1 aseptic field, 1 dressing tray, 1 paper towel, 2 non-woven swabs, 5 non-woven balls, 1 yellow forcep, 2 blue forceps	160
08-666NP	1 aseptic field, 1 dressing tray, 1 paper towel, 2 gauze swabs, 5 cotton balls, 1 yellow forcep, 2 blue forceps	160
07-889P	1 aseptic field, 1 dressing tray, 1 paper towel, 3 gauze swabs, 6 cotton balls, 1 yellow forcep, 2 blue forceps	120
06-782	1 aseptic field, 1 dressing tray, 4 non-woven swabs, 2 blue forceps	160

Tray 2: Dressing Packs Peel Packs		
RE-ORDER	DESCRIPTION	PER CARTON
08-118P	1 aseptic field, 1 dressing tray, 1 laminated towel, 6 non-woven swabs, 1 yellow forcep, 2 blue forceps	160
08-333P	1 aseptic field, 1 dressing tray, 1 paper towel, 6 non-woven swabs, 1 yellow forcep, 2 blue forceps	160
08-709	1 aseptic field, 1 dressing tray, 1 laminated towel, 6 gauze swabs, 1 yellow forcep, 2 blue forceps	160
07-190	1 aseptic field, 1 dressing tray, 1 disposable bag white, 1 paper towel, 2 non-woven swabs, 1 yellow forcep, 2 blue forceps, 5 non-woven balls	160
08-625	1 aseptic field, 1 dressing tray, 1 laminated towel, 2 gauze swabs, 6 non-woven swabs, 1 yellow forcep, 2 blue forceps	170

For more information please contact your Business Development Manager:

(02) 8718 2888 or visit multigate.com.au