

Eisner™

A division of Multigate
Medical Products Pty Ltd

Enhance your every move







Eisner™ laparoscopic
instruments deliver
unmatched durability,
comfort and efficiency
to support a seamless
performance every time.



POWERING
LAPAROSCOPIC
SURGERY
SINCE 1993



CONTENTS

	Handles	6		Needleholders	40
	'C' Series	8		Handles	41
	'G' Series	9		Shafts	42
	'R' Series	10		Inserts	42
	'S' Series	11		Accessories	43
	Inserts - Forceps (5mm)	12		Other Instruments	44
	Biopsy	14		Bulldog Clamps	45
	Ovary/Tubal	14		Retractors	46
	Clamps	15		Cholangiography	46
	Curved Vascular Clamps	16		Knot Guide	47
	Dissectors	17		Knot Pushers	47
	Fenestrated Dissectors	18		Myoma Drills	48
	Extractors	20		Palpation Probes	48
	Graspers	21		Puncture Needles	49
	Grasper Scissors	23		Accessories & Maintenance	50
	Inserts - Forceps (10mm)	24		Silicon Valves	51
	Biopsy	24		Sealing Caps	52
	Clamps	24		Cleaning Brush	52
	Dissectors	24		Lubricating Oil	52
	Extractors	25		Single Use Instruments	54
	Graspers	25		Scissors	56
	Inserts - Scissors	26		Scissor Insert & Sheath	56
	Scissors with Insulated Tube	28		Forceps	57
	Scissors - Curved Blade	28		Monopolar Diathermy	58
	Scissors - Straight Blade	29		Suction Irrigation Devices	58
	Other Scissors	29		Specimen Retrieval	59
	Insulated Tubes, Sheaths & Locking Nuts	30			
	Safety-Lip	31			
	Ezi-Clean	31			
	Locking Nuts 5mm	31			
	Molded Hardcoat	32			
	Locking Nuts 10mm	32			
	Diathermy - Monopolar	34			
	Electrodes	35			
	Electrodes Ceramatip®	35			
	Cables	35			
	Diathermy - Bipolar	36			
	Handles	37			
	Shafts/Tubes	37			
	Non-stick Forceps	37			
	Electrodes	38			
	Cables	38			




DESIGN

Our instruments are designed for optimal balance, feel, comfort and responsiveness to deliver a seamless connection between the hand and the surgical site.

SAFETY & RELIABILITY

The craftsmanship and the quality of materials we select combine to deliver high performance, electro-safe instruments that deliver unrivalled durability, longevity and reliability.



- **Easy to re-process** Our instruments are designed to simplify re-processing.

- **Clinical performance** Our hand-held instruments are modular allowing you to build cost effective instrument trays whilst maximising patient safety and clinical performance.



Safety Additional insulation improves protection for the user.

Comfort Rounded — ergonomic handle improves comfort and reduces fatigue.

Precision Ratchet mechanism provides excellent, durable tissue and grip control.

Eisner™ | eisner@multigate.com.au

Handles - 'C' Series

Our composite 'C' series handles come as standard, ratchet, multi-function and swing away ratchet versions.

'C' Series Standard Composite Handle		
RE-ORDER	FEATURES	MIN QTY
CFS	Standard (Non-Ratchet)	1



'C' Series Ratchet Composite Handle		
RE-ORDER	FEATURES	MIN QTY
CFR	Ratchet	1



'C' Series Multifunction Ratchet Composite Handle		
RE-ORDER	FEATURES	MIN QTY
CFM	Multifunction Standard/Ratchet	1



'C' Series Swing Away Ratchet Composite Handle		
RE-ORDER	FEATURES	MIN QTY
CFZ	Swing Away Ratchet	1



Handles - 'G' Series

Our composite 'G' series handles come as standard, ratchet, and ratchet with on/off functionality.

'G' Series Standard Composite Handle		
RE-ORDER	FEATURES	MIN QTY
GSS	Standard (Non-Ratchet)	1



'G' Series Ratchet Composite Handle		
RE-ORDER	FEATURES	MIN QTY
GSR	Ratchet (Push Button)	1



'G' Series On/Off Ratchet Composite Handle		
RE-ORDER	FEATURES	MIN QTY
GSA	Convertible Standard/Ratchet	1



Handles - 'R' Series

Our composite 'R' series handles come in standard and ratchet versions.

'R' Series Standard Composite Handle		
RE-ORDER	FEATURES	MIN QTY
RNR	Standard (Non-Ratchet)	1



'R' Series Ratchet Composite Handle		
RE-ORDER	FEATURES	MIN QTY
RRH	Ratchet	1



Handles - 'S' Series

Our stainless steel 'S' series handles come in standard, ratchet and spring return versions.

'S' Series Standard Steel Handle		
RE-ORDER	FEATURES	MIN QTY
SSS	Standard (Non-Ratchet)	1



'S' Series Ratchet Steel Handle		
RE-ORDER	FEATURES	MIN QTY
SSR	Ratchet	1



'S' Series Spring Return Steel Handle		
RE-ORDER	FEATURES	MIN QTY
SRH	Spring Return Ratchet	1



'S' Series Swing Away Ratchet Steel Handle		
RE-ORDER	FEATURES	MIN QTY
SAR	Swing Away Ratchet	1





INSERTS - FORCEPS

Eisner™ inserts are tough performers designed to deliver high fidelity control without compromising on their durability, lifespan and reliability.

- **Easy to manage** Our inserts are designed to self-open for sterilisation.

- **Durable** No weld design ensures long-term performance.








Working Length Inserts available in standard and bariatric lengths.

Performance Broad range of precision jaws to support all types of procedure.

Inserts - Forceps (5mm)

We offer a broad range of jaw types in both standard and bariatric lengths to cover all minimally invasive procedures.

All our inserts are fully compatible and inter-changeable and designed for durable, reliable performance.

Biopsy	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Biopsy Forceps	400082	401082	11mm Jaw	Double Action Sharp Oval Cups	1
	5mm Biopsy Pin Forceps	400521	401521	10mm Jaw	Single Action Sharp Oval Cups 2 Pins	1
	5mm Biopsy Punch Insert	400100	401100	11mm Jaw	Single Action Thru-Cutting Mechanism	1
	5mm Blakesley Biopsy Forceps	400086	401086	10mm Jaw	Single Action Sharp Elongated Fenestrated Oval Cups	1
	5mm Minibite Biopsy Forceps	400080	401080	8mm Jaw	Double Action Sharp Short Oval Cups 2 Pins	1







Ovary/Tubal	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Blunt Nose Tubal Forceps	400046	401046	20mm Jaw	Double Action Ovary/Tubal Forceps Non Serrated Inner Jaw	1
	5mm Tubal Forceps	400045	401045	16mm Jaw	Double Action Fine Horizontal Serrations	1

INSERTS - FORCEPS (5MM)









Clamps	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Debakey Forceps	400102	401102	20mm Jaw	Double Action 2 into 3 Debakey Serrations Slim Line Profile	1
	5mm Debakey Forceps	400103	401103	30mm Jaw	Double Action 2 into 3 Debakey Serrations Slim Line Profile	1
	5mm Debakey Forceps	400104	401104	40mm Jaw	Double Action 2 into 3 Debakey Serrations Slim Line Profile	1
	5mm Dorsey Fenestrated Forceps	400287	401287	45mm Jaw	Double Action Heavy Horizontal Serrations Long Single Fenestration	1
	5mm Hunter Jaeger Forceps	400012	401012	21mm Jaw	Double Action Debakey Serrations Long Single Fenestration	1
	5mm Marques DA Double Fen. Forceps	400311	401311	25mm Jaw	Double Action Medium Horizontal Serrations Double Fenestration	1
	5mm Marques DA Triple Fen. Forceps	400313	401313	35mm Jaw	Double Action Medium Horizontal Serrations Triple Fenestration	1
	5mm Marques SA Double Fen. Forceps	400310	401310	25mm Jaw	Single Action Medium Horizontal Serrations Double Fenestration	1

INSERTS - FORCEPS (5MM)






Clamps (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Marques SA Triple Fen. Forceps	400312	401312	35mm Jaw	Single Action Medium Horizontal Serrations Triple Fenestration	1
	5mm Marques Upcurve Forceps	400314	401314	25mm Jaw	Single Action Medium Horizontal Serrations Double Fenestration	1

Curved Vascular Clamps	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Crafoord 35 Clamp Forceps	400204	401204	35mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	5mm Crafoord 50 Clamp Forceps	400206	401206	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	5mm Debakey Curved 50 Clamp Forceps	400220	401220	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	5mm Debakey S-shape 50 Clamp Forceps	400221	401221	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	5mm Satinsky 40 Clamp Forceps	400202	401202	40mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1
	5mm Satinsky 50 Clamp Forceps	400200	401200	50mm Jaw	Single Action Atraumatic 1 into 2 Debakey Serrations	1

INSERTS - FORCEPS (5MM)








Dissectors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Dolphin Cup Forceps	400015	401015	15mm Jaw	Double Action Medium Horizontal Serrations with Cup	1
	5mm Dolphin Forceps	400021	401021	14mm Jaw	Double Action Medium Horizontal Serrations	1
	5mm Duckbill Rounded Forceps	400087	401087	10mm Jaw	Double Action Parallel/Aligned Jaw Closure Medium Horizontal Serrations	1
	5mm Geegee Forceps	400030	401030	38mm Jaw	Single Action Medium Horizontal Serrations	1
	5mm Kelly Baby Criss Cross Forceps	400733	401733	15mm Jaw	Double Action Criss Cross Serrations	1
	5mm Kelly Criss Cross Forceps	400732	401732	20mm Jaw	Double Action Criss Cross Serrations	1
	5mm Maryland Debakey Forceps	551712	—	29mm Jaw	Double Action 1 into 2 Debakey Serrations Fine Teeth at Tip	1
	5mm Maryland Long Forceps	400660	401660	24mm Jaw	Double Action Criss Cross Serrations	1

INSERTS - FORCEPS (5MM)









Dissectors (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Maryland Standard Forceps	400731	401731	20mm Jaw	Double Action Fine Horizontal Serrations	1
	5mm Maryland Grasper Dissector Forceps	400661	401661	17mm Jaw	Double Action Deep Horizontal Serrations	1
	5mm Mixer Forceps	400735	401735	13mm Jaw	Double Action Parallel/Aligned Jaw Closure Fine Criss Cross Serrations	1
	5mm Pointed Forceps	400016	401016	22mm Jaw	Double Action Criss Cross Serrations	1
	5mm Reddick Olsen Forceps	400342	401342	13mm Jaw	Double Action Medium Horizontal Serrations	1

Fenestrated Dissectors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Bullnose Fenestrated Forceps	400040	401040	10mm Jaw	Double Action Medium Horizontal Serrations Single Fenestration	1
	5mm Fenestrated Dissector Forceps	400088	401088	13mm Jaw	Double Action Fine Horizontal Serrations Single Elongated Fenestration	1









INSERTS - FORCEPS (5MM)

Fenestrated Dissectors (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Fenestrated Strong Forceps	400132	401132	8mm Jaw	Double Action Medium Horizontal Serrations Single Elongated Fenestration	1
	5mm Johans Fenestrated Forceps	400285	401285	17mm Jaw	Double Action Medium Horizontal Serrations Spacious Fenestration	1
	5mm Johans Long DA Fen. Forceps	400291	401291	40mm Jaw	Double Action Super Fine Horizontal Serrations Long Spacious Fenestration	1
	5mm Johans Long SA Fen. Forceps	400288	401288	40mm Jaw	Single Action Super Fine Horizontal Serrations Long Spacious Fenestration	1
	5mm Johans Loop DA Forceps	400289	401289	21mm Jaw	Double Action Medium Horizontal Serrations Spacious Fenestration	1
	5mm Johans Loop SA Forceps	400284	401284	20mm Jaw	Single Action Medium Horizontal Serrations Spacious Fenestration	1
	5mm Wave Fenestrated Forceps	400052	401052	14mm Jaw	Double Action Fine Horizontal Serrations Single Long Fenestration	1









INSERTS - FORCEPS (5MM)

Extractors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Cats Paw Tenaculum Forceps	400048	401048	20mm Jaw	Double Action 2 into 3 End Teeth 2 Pair Outer Side Hook Teeth	1
	5mm Davis Claw (Serrated Tooth) Forceps	400717	401717	22mm Jaw	Double Action 2 Into 3 End Teeth Deep Horizontal Grooves	1
	5mm Eistige Forceps	400013	401013	22mm Jaw	Double Action Deep Cup Rippled Edges	1
	5mm Eistige SA Forceps	400033	401033	22mm Jaw	Single Action Deep Cup Rippled Edges	1
	5mm Gripper Forceps	400017	401017	21mm Jaw	Double Action Coarse Serrations Long Fenestration	1
	5mm Interlocker Forceps	400217	401217	22mm Jaw	Double Action Interlocking Smooth Raised Teeth Long Fenestration	1
	5mm Ripple Claw (St Andrews) Forceps	400715	401715	23mm Jaw	Double Action 2 into 4 Aggressive End Teeth	1
	5mm Tenaculum Forceps	400364	401364	22mm Jaw	Double Action Super Fine Claw Tip Smooth Inner Jaw	1








INSERTS - FORCEPS (5MM)

Graspers	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Alligator Retrograde Forceps	400010	401010	15mm Jaw	Double Action Strong Retrograde Serrations	1
	5mm Allis 35 Forceps	400197	401197	35mm Jaw	Double Action 4mm Wide Double Row Platform Smooth Inner Jaw	1
	5mm Allis Forceps	400089	401089	13mm Jaw	Double Action Serrated End Platform Smooth Inner Jaw	1
	5mm Allis SA Forceps	400011	401011	13mm Jaw	Single Action Serrated End Platform Smooth Inner Jaw	1
	5mm Allis Long Forceps	400330	401330	28mm Jaw	Double Action Fine Toothed Tip Smooth Inner Jaw	1
	5mm Babcock Duval Forceps	400130	401130	25mm Jaw	Double Action Horizontal Serrations on Medium Platform Elongated Fenestration	1
	5mm Babcock SA Forceps	400251	401251	12mm Jaw	Single Action Fine Horizontal Serrations Long Fenestration	1
	5mm Babcock Super Fine Paddle Forceps	400250	401250	30mm Jaw	Double Action 5mm Platform with Horizontal Serrations Smooth Inner Jaw	1

INSERTS - FORCEPS (5MM)


Graspers (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Babcock Traditional Forceps	400252	401252	20mm Jaw	Double Action Distal End Fenestration Small Platform with Horizontal Serrations	1
	5mm Bonnies Forceps	400009	401009	15mm Jaw	Double Action 1 into 2 Teeth Fine Horizontal Serrations	1
	5mm Bullnose Flat Face Forceps	400355	401355	13mm Jaw	Double Action Medium Horizontal Serrations	1
	5mm Bullnose Forceps	400127	401127	13mm Jaw	Double Action Medium Horizontal Serrations	1
	5mm Bullnose Long Forceps	400002	401002	22mm Jaw	Double Action Medium Horizontal Serrations	1
	5mm Bullnose Short Jaw Forceps	400001	401001	11mm Jaw	Double Action Medium Horizontal Serrations	1
	5mm Crocodile Forceps	400019	401019	17mm Jaw	Double Action Precision Pyramidal Teeth Small Triple Fenestration	1
	5mm Kocher Forceps	400004	401004	15mm Jaw	Double Action 2 into 3 End Teeth Fine Horizontal Serrations	1


INSERTS - FORCEPS (5MM)






Graspers (continued)	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm McKernan Medium Forceps	400025	401025	20mm Jaw	Double Action Traditional McKernan Style	1
	5mm McKernan Standard Forceps	400136	401136	15mm Jaw	Double Action Traditional McKernan Style	1
	5mm McKernan Long Forceps	400137	401137	22mm Jaw	Double Action McKernan Style Long Single Fenestration	1
	5mm Rat Tooth Forceps	400036	401036	12mm Jaw	Double Action 2 into 4 Teeth Medium Horizontal Serrations	1
	5mm Pin Grasper Forceps	400042	401042	12mm Jaw	Double Action 1 into 2 End Teeth Smooth Inner Jaw	1
	5mm Tungsten Carbide (TC) Grasper Forceps	400029	401029	16mm Jaw	Double Action Tungsten Carbide Jaw Pyramidal Serrations	1
	5mm Cobra Forceps	400201	401201	13mm Jaw	Double Action 2 into 4 Sharp End Teeth	1
Grasper Scissors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Knot Pusher Scissor Insert	400331	401331	16mm Jaw	Double Action Knot Pusher Tip and Scissors in Unique Molding	1

Inserts - Forceps (10mm)


The density of our 10mm insulation tubes is tripled to minimize splitting, cracking and deformity and their single piece design with no welded joints delivers heavy duty, reliable performance every time.


Biopsy	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	10mm Spoon Forceps	400720	401720	27mm Jaw	Single Action Sharp Large Cups Rounded End Profile	1






Clamps	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	10mm Nussbaum Intestinal Forceps	400781	401781	90mm Jaw	Double Action Full Length Vertical Serrations Additional Fine Cross-Serrations on Leading Section	1

Dissectors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	10mm Geegee Forceps	400076	401076	52mm Jaw	Single Action Medium Horizontal Serrations	1
	10mm Mixer 90° Forceps	400739	401739	28mm Jaw	Double Action Medium Horizontal Serrations	1
	10mm Overholt Criss Cross Forceps	400741	401741	25mm Jaw	Double Action Medium Criss Cross Serrations	1
	10mm Overholt-Mixer Stepdown Forceps	400737	401737	30mm Jaw	Double Action Deep Horizontal Serrations	1
	10mm Mixer Right Angle Forceps	400736	401736	24mm Jaw	Double Action Fine Criss Cross Serrations	1

INSERTS - FORCEPS (10MM)

Extractors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	10mm Claw (Mother in-law) Forceps	400710	401710	30mm Jaw	Single Action Aggressive 2 into 3 End Teeth Smooth Inner Jaw	1


	10mm Tenaculum Forceps	400753	401753	35mm Jaw	Double Action Fine Claw Tip Smooth Inner Jaw	1
---	------------------------	--------	--------	----------	--	---

Graspers	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	10mm Anvil Grasping Forceps	400295	401295	30mm Jaw	Double Action Specially Shaped to Hold the Anvil of a Circular Stapler	1
	10mm Babcock Forceps	400078	401078	30mm Jaw	Double Action Medium Platform with Fine Vertical Serrations V-shaped Fenestration	1
	10mm Collin-Duval Grasping Forceps	400084	401084	32mm Jaw	Double Action Criss Cross Serrations on Platform Long Broad Fenestration	1
	10mm Duval Long Forceps	400280	401280	50mm Jaw	Double Action Deep Vertical Serrations on Platform Short V-shaped Fenestration	1
	10mm Pennington-Collin Forceps	400085	401085	38mm Jaw	Double Action Deep Horizontal Serrations Spacious V-shaped Fenestration	1





INSERTS - SCISSORS

We offer a range of straight and curved blade scissor inserts.






Inserts - Scissors (5mm)



Scissors with Insulated Tube	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Metz Curved Scissor (Tri-Cut) with Insulated Tube & L'nut	300-600-177	300-600-177.45	15mm Blade	Double Action Medium Curve Blade Sheath and Locking Nut Included	1

Curved Blade	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Curved Muticut Scissors	MC2	—	11mm Blade	Double Action Mild Curve Maximum of 10 Procedure Use Recommended	1
	5mm Curved Multicut Scissors (Pack of 5)	MC2000	—	11mm Blade	Double Action Mild Curve Blade	1 Pack

	5mm Metz Strong Curve Baby Scissors Insert	552030	—	14mm Blade	Double Action Strong Curve Metzenbaum Blade	1
	5mm Metz Strong Curve Scissors Insert	552112	—	16mm Blade	Double Action Strong Curve Metzenbaum Blade	1
	5mm Metz Medium Curve Scissor Insert	400067	401067	12mm Blade	Double Action Medium Curve Metzenbaum Blade	1
	5mm Metz. Strong Curve NPHM Scissor Insert	400072	401072	15mm Blade	Double Action Non-protruding Hinge Mechanism	1
	5mm Micro Curved Scissors Insert	400572	401572	7mm Blade	Double Action Short Petite Curved Blade	1

INSERTS - SCISSORS (5MM)


Straight Blade	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Metz Straight Scissor Insert	400069	401069	12mm Blade	Double Action Straight Metzenbaum Blade	1
	5mm Micro Straight Scissors Insert	400570	401570	7mm Blade	Single Action Short Blade Wide Opening Feature	1
	5mm Metz Straight Scissors Insert	400068	401068	15mm Blade	Double Action Straight Longer Metzenbaum Blade	1
	5mm Straight Multicut Scissors	MC1	—	11mm Blade	Double Action Straight Blade	1
	5mm Straight Multicut Scissors (Pack of 5)	MC1000	—	11mm Blade	Double Action Straight Blade	1 Pack

Other Scissors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	5mm Hook Scissor Insert	400060	401060	8mm Blade	Single Action Hook Style Blade	1
	5mm Knot Pusher Scissor Insert	400331	401331	16mm Blade	Double Action Knot Pusher Tip & Scissor in Unique Molding	1




INSULATED TUBES, SHEATHS & LOCKING NUTS

Our insulated tubes are designed for easy re-processing and to reduce the potential for fluids, debris and bio-burdens to enter under the insulation.



10MM INSULATED TUBE

Durable Insulation material density increased three-fold to minimise splitting, cracking or deformity.



5MM INSULATED TUBE

Ease Our Ezi-Clean™ Insulated Tubes have a removable sheath for easy cleaning and simple assembly.



Our range of coloured locking nuts help you to code your instruments trays to specific laparoscopic procedures.


Insulated Tubes, Sheaths & Locking Nuts (5mm)






Our insulated tubes are made from surgical grade stainless steel and have a Polyphenylsulfone coating designed to resist heat and chemicals.

Tubes and sheaths are available in standard and bariatric lengths and in both 5mm and 10mm diameters.

Safety-Lip	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	FEATURES	MIN QTY
	5mm Safety-Lip Insulated Tube with locking nut	IT-05/33	IT-05/45	Steel Safety-lip Sealed Tip	1






Ezi-Clean	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	FEATURES	MIN QTY
	5mm Ezi-Clean Insulated Tube with locking nut	IT-05/33.EZC	IT-05/45.EZC	Insulation Tube is Removable Re-usable & Replaceable	1

	Black Insulation Sheath Only (Pack of 5)	EZC-33	EZC-45	Replaceable Insulation Sheaths	1 Pack
---	--	--------	--------	--------------------------------	--------

Locking Nuts 5mm	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	5mm Eisner Locking Nut - Standard Black	ELN-05SB	Ø5mm	Standard Black	1
	5mm Eisner Locking Nut - Blue	ELN-05B	Ø5mm	Blue	1
	5mm Eisner Locking Nut - Green	ELN-05G	Ø5mm	Green	1
	5mm Eisner Locking Nut - Red	ELN-05R	Ø5mm	Red	1
	5mm Eisner Locking Nut - Yellow	ELN-05Y	Ø5mm	Yellow	1

Insulated Tubes, Sheaths & Locking Nuts (10mm)

Molded Hardcoat	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	FEATURES	MIN QTY
	10mm Molded Hardcoat Insulation Tube with locking nut	IT-10/33	IT-10/45	Thickened Hard-Coat Insulation	1

Locking Nuts 10mm	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	10mm Eisner Locking Nut - Standard Black	ELN-10SB	Ø10mm	Standard Black	1
	10mm Eisner Locking Nut - Blue	ELN-10B	Ø10mm	Blue	1
	10mm Eisner Locking Nut - Green	ELN-10G	Ø10mm	Green	1
	10mm Eisner Locking Nut - Red	ELN-10R	Ø10mm	Red	1
	10mm Eisner Locking Nut - Yellow	ELN-10Y	Ø10mm	Yellow	1

We work at the intersection between the surgeon and their patient, designing high performance, responsive laparoscopic instruments that offer a seamless connection between the hand and the surgical site.




DIATHERMY - MONOPOLAR

Eisner™ Monopolar Electrodes are uncompromising on patient and clinician safety. Every instrument must meet exacting standards that far exceed the Australian standard for electro-safety.



Safety Overlapping electro-safe cable connection.





Durable Our light-weight Ceramatip™ electrodes feature a hardened durable crystalline ceramic tip.

Diathermy - Monopolar

Eisner™ Electrodes are uncompromising on patient and clinician safety. Every instrument we make must meet our own exacting standards which far exceed Australian standards. Our 300-100 series electrodes are made from a hard coat PPSU polymer and come with an Electro-Safe grip & plug on their distal end.

Our light-weight Ceramatip™ 'L'-Hook and Spatula feature a hardened durable crystalline-ceramic tip, a stop-roll handle moulding, quality hard PEEK shaft insulation and an overlapping electro-safe cable connection.

Electrodes	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	FEATURES	MIN QTY
	Spatula Electrode	300-100-131	300-101-131	Strong Curved Spade	1
	J-Hook Electrode	300-100-132	—	Smooth L-Hook Tip	1
	L-Hook Electrode	300-100-133	300-101-133	Premium Contoured L-Hook	1
	Button Electrode	300-100-134	—	Smooth Cylindrical Contour	1

Electrodes Ceramatip®	DESCRIPTION	RE-ORDER 34cm	RE-ORDER 45cm	FEATURES	MIN QTY
	Ceramatip® L-Hook Electrode with Ceramic Tip	300-163-420	300-163-620	Standard L-Hook Design	1
	Ceramatip® Spatula Electrode with Ceramic Tip	300-163-430	300-163-630	Long Slender Curved Spade	1

Cables	DESCRIPTION	RE-ORDER 3m	RE-ORDER 5m	SIZE	FEATURES	MIN QTY
	3m & 5m Monopolar Diathermy Cable	700-360-031	700-360-051	3m 5m	4mm Round Plug Valleylab Compatible	1



DIATHERMY - BIPOLAR

Eisner™ offers a range of high-quality Bipolar electrodes for use with our Powergrip contour handle. Our cables meet the highest safety standards and come with a twin moulded plug compatible with Valleylab, EMC & Lamidey.




Diathermy - Bipolar






Bipolar Handle		
RE-ORDER	FEATURES	MIN QTY
300-860-134	Powergrip contour	1



Shafts/Tubes		DESCRIPTION	RE-ORDER 34cm	RE-ORDER 45cm	SIZE 34cm	SIZE 45cm	FEATURES	MIN QTY
		Bipolar Tube/Shaft	300-860-141	300-860-147	34cm WL	45cm WL	Quality Surgical Grade Stainless Steel	1

Non-Stick Forceps		DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
		Bipolar Non-Stick Forceps	700-082-200	200mm	Bayonet Design Straight 2mm Wide Tip	1

DIATHERMY - BIPOLAR

Electrodes	DESCRIPTION	RE-ORDER 34cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	Bipolar Slightly Curved Grasper Electrode	300-860-011	300-860-021	20mm Jaw	Double Action Mild Curve Electrode Insulated to Tip of Grasper	1
	Bipolar Delta Fenestrated Electrode	300-860-014	300-860-024	13mm Jaw	Double Action Single Fenestration	1
	Bipolar Curved Grasper Electrode	300-860-015	300-860-025	18mm Jaw	Double Action Mild Curve Fine Horizontal Serrations Small Fenestration	1
	Bipolar TC-Metzenbaum Scissors Electrode	300-860-018	300-860-028	15mm Jaw	Double Action Tungsten Carbide Mild Curve Blade	1
	Bipolar Maryland Electrode	300-860-019	300-860-029	21mm Jaw	Double Action Mild Slender Curve Fine Horizontal Serrations	1

Cables	DESCRIPTION	RE-ORDER 3m	RE-ORDER 5m	FEATURES	MIN QTY
	3m & 5m Bipolar Diathermy Cable	700-351-031	700-351-051	Round Mold Flat Plug Connector	1
	3m & 5m Bipolar Diathermy Cable	700-351-131	700-351-151	Square Mold Flat Plug Connector	1
	3m & 5m Bipolar Diathermy Cable	700-351-231	700-351-251	2 Pin Plug Angled Molding Valleylab Compatible	1

We have always felt the highest duty of care to you and your patients. With every iteration of our products we aim to set the standard in safety and performance.



NEEDLEHOLDERS

Our Genesis Needleholder series offers ergonomic axial and pistol style handles together with a range of jaw inserts as part of a modular three-piece dismantlable design that simplifies re-processing whilst delivering optimal comfort, tissue grip and needle guidance during suturing.

Flushable
cleaning
ports

- Three-piece dismantlable design simplifies re-processing
- Ergonomic design
- Available in 33cm and 45cm lengths
- Range of different jaw configurations

Available in
stainless steel
and anodised
aluminium

Needleholders - Handles

Genesis Handle I - Stainless Steel Axial Style

RE-ORDER	FEATURES	MIN QTY
300-480-002	Axial Style	1



Genesis Handle II - Lightweight Pistol Style

RE-ORDER	FEATURES	MIN QTY
300-480-001	Lightweight Pistol Style	1



Genesis Handle III - Lightweight Axial Style

RE-ORDER	FEATURES	MIN QTY
300-480-000	Lightweight Axial Style	1





Needleholders - Shafts



Shafts	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	330mm Genesis Shaft/Tube	300-480-033	330mm WL	Stainless Steel Tube with Stainless Steel Locking Nut	1
	450mm Genesis Shaft/Tube	300-480-045	450mm WL	Stainless Steel Tube with Stainless Steel Locking Nut	1

Needleholders - Inserts






Inserts	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	Genesis Insert Straight Jaw	300-405-130	300-406-130	12mm Jaw	Straight Tapered Single Action Jaw Fine Criss Cross Serration	1
	Genesis Insert Left Curve Jaw	300-405-140	300-406-140	12mm Jaw	Left Curve Single Action Jaw Fine Criss Cross Serrations	1
	Genesis Insert Right Curve Jaw	300-405-141	300-406-141	12mm Jaw	Right Curve Single Action Jaw Fine Criss Cross Serrations	1
	Genesis Insert Diamond Dusted - Self Adjusting Jaw	300-405-142	300-406-142	12mm Jaw	Cupped Lower Single Action Jaw with Self-adjusting Upper Jaw	1



Needleholders - Accessories



Accessories	DESCRIPTION	RE-ORDER	FEATURES	MIN QTY
	Genesis Steel Locking nut	DR30-50-1055S	Stainless Steel	1
	Flush Port/ Luer lock seal, rubber	300-400-901	Black Rubber	1
	Flush Port/ Luer lock seal, stainless steel	300-999-010	Stainless Steel	1





OTHER
INSTRUMENTS



Other Instruments - Bulldog Clamps

Bulldog Clamp - Klein	DESCRIPTION	RE-ORDER 33cm	SIZE	FEATURES	MIN QTY
	Endoscopic Klein® Bulldog Clamp Applier	300-410-775	Ø10mm	Applier for Detachable Clamps	1
	Endoscopic Klein® Bulldog Bowel Straight Clamp	300-100-771	45mm Jaw	Straight Mid-length Clamp	1
	Endoscopic Klein® Bulldog Bowel Straight Clamp	300-100-772	60mm Jaw	Straight Long-length Clamp	1
	Endoscopic Klein® Bulldog Vascular Straight Clamp	300-100-773	17mm Jaw	Straight Short-length Clamp	1
	Endoscopic Klein® Bulldog Bowel Straight Clamp	300-100-774	25mm Jaw	Straight Mid-length Clamp	1
	Endoscopic Klein® Bulldog Vascular Curved Clamp	300-100-776	17mm Jaw	Mild Curve Short-length Clamp	1
	Endoscopic Klein® Bulldog Vascular Curved Clamp	300-100-777	25mm Jaw	Mild Curve Mid-length Clamp	1
	Endoscopic Klein® Bulldog Bowel Curved Clamp	300-100-778	45mm Jaw	Mild Curve Long-length Clamp	1
	Endoscopic Klein® Bulldog Bowel Curved Clamp	300-100-779	60mm Jaw	Mild Curve Extra Long-length Clamp	1
	Endoscopic Klein® Bulldog Satinsky Back Curved Trocar Clamp	300-100-786	60mm Jaw	Long Back Curved Trocar Clamp	1
	Endoscopic Klein® Bulldog Satinsky Forward Curved Trocar Clamp	300-100-789	60mm Jaw	Long Forward Curved Trocar Clamp	1



Other Instruments - Retractors

Retractor	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	10mm Fan Retractor Dismantlable Elevator	300-065-100	380mm WL	5 Finger Fan Spread Articulating Head Dismantlable	1
	S-shaped Tissue Retractor	300-098-105	180mm	Stainless Steel S-shaped Tissue Retractor Set of 2 pcs	1



Other Instruments - Cholangiography



Cholangiography	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	Rubber Seal/Silicone Bung	301-701-100		Pale Grey Port Bung Specifically for Cholangiogram Forceps	1
	5mm Cholangiogram Forceps 330mm W/L	300-100-760	8mm Jaw	For 6 Fr Gauge Catheter (2mm)	1
					



Other Instruments - Knot Guide

Knot Guide	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	FEATURES	MIN QTY
	Knot Guide	300-030-000	300-030-000.45	Long Cylindrical Form with Slit	1




Other Instruments - Knot Pushers

Knot Pushers	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	FEATURES	MIN QTY
	Knot Pusher - Open Face	300-030-001	300-030-003	Ø4mm Classic Knot Pusher with Slit and Eye	1
	Knot Pusher - Closed face	300-030-004	300-030-006	Ø4mm Classic Knot Pusher with Slit and Eye	









Other Instruments - Myoma Drills


Myoma Drills	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	5mm Myoma Drill/Screw	300-032-015	330mm WL	Stainless Steel 1pce Instrument Cork Screw End	1
	5mm Myoma Drill/Screw	300-042-005	420mm WL	Stainless Steel 1pce Instrument Cork Screw End	1
	10mm Myoma Drill/Screw	300-033-010	330mm WL	Stainless Steel 1pce Instrument Cork Screw End	1
	11mm Myoma Drill/Screw	300-033-011	330mm WL	Stainless Steel 1pce Instrument Cork Screw End	1

+ Other Instruments - Palpation Probes

Palpation Probes	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	5mm Palpation Probe with Markings	300-193-305	300mm WL	With Graduated Depth Markings and Small Handles	1
	5mm Palpation Probe with Markings	300-193-335	330mm WL	With Graduated Depth Markings and Small Handles	1
	5mm Palpation Probe with Markings	300-193-405	400mm WL	With Graduated Depth Markings and Small Handles	1
	5mm Palpation Probe with Markings	300-193-455	450mm WL	With Graduated Depth Markings and Small Handles	1
	10mm Palpation Probe with Markings	300-193-410	400mm WL	With Graduated Depth Markings and Small Handles	1



Other Instruments - Puncture Needles


Puncture Needles	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	Puncture Needle with Luer - Lock Adaptor	300-194-307	Ø0.7mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-309	Ø0.9mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-312	Ø1.2mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-314	Ø1.4mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-318	Ø1.8mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-412	Ø1.2mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-414	Ø1.4mm Tip	Ø5mm - 20mm Needle Length	1
	Puncture Needle with Luer - Lock Adaptor	300-194-418	Ø1.8mm Tip	Ø5mm - 20mm Needle Length	1





ACCESSORIES & MAINTENANCE



We offer a range of accessories and maintenance equipment to help you maintain optimal tray performance.



Accessories & Maintenance - Sealing Caps

Sealing Caps	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	10mm Sealing Cap Green	300-025-710	Ø10.5mm ID	For use with 10mm Hasson Cone or 12.5mm Port	1
	11mm Sealing Cap Grey	300-025-711	Ø11mm ID	For use with 11mm Hasson Cone	1
	10mm Sealing Cap Purple	300-025-910	Ø8mm ID	For use with 10mm Port	1
	11mm Sealing Cap Blue	300-025-911	Ø8.5mm ID	For use with 11mm Port	1
	5-5.5mm Sealing Cap Green	300-025-955	Ø3.5mm ID	For use with 5-5.5mm Port	1
	Spring Cap for all size Ergolite Ports	300-025-999		For use with all Size Ergolite Ports Screw Thread	1
	Rubber seal/ Silicone Bung	301-701-100		Pale Grey Specifically for Cholangiogram Forceps	1
	Flush Port/ Luer Lock Seal, Rubber	300-400-901		Black Rubber	1
	Flush Port/ Luer Lock Seal, Stainless Steel	300-999-010		Stainless Steel	1


Accessories & Maintenance - Silicon Valves

Silicone Valves	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	Silicone Valve for Ø3-5.5mm Ports/Trocars (Pack of 5)	300-025-355		Fits Ø3mm - 5.5mm Port	1 Pack
	Silicone Valve for Ø10-12.5mm Ports/Trocars (Pack of 5)	300-025-111		Fits Ø10mm - 12.5mm Port	1 Pack

Accessories & Maintenance - Cleaning Brush

Cleaning Brush	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	Lap Cleaning Brush for 5mm Instrument	300-135-905	470mm WL	Black Ring Handle Standard Brush for 5mm Instrument	1
	Lap Cleaning Brush for 10mm Instrument	300-135-910	450mm WL	Black Ring Handle Standard Brush for 10mm Instrument	1

Accessories & Maintenance - Lubricating Oil

Lubricating Oil	DESCRIPTION	RE-ORDER	SIZE	FEATURES	MIN QTY
	Instrument Oil 50ml Bottle	700-900-160	50ml	Squeezable Bottle (not Aerosol)	1

Instrument Care

Contact us to book a **FREE** assessment of your trays.

Call **(03) 9728 4770**
or visit eisner.com.au/services



2

SINGLE USE INSTRUMENTS

We understand the impact inferior equipment can have on your performance, which is why we have worked closely with surgeons in the development of our Cadence™ single use laparoscopic instruments.

From the intuitive, versatile and ergonomic design of our Suction Irrigator, to the unique one-handed deployment of our Specimen Retrieval Device, Cadence™ offers a single use solution with a wide

range of superior features designed to deliver the precision, comfort and control that you expect from reusable laparoscopic instruments.

Performance Precision tip provides constant, precise irrigation flow without compromising suction.



Intuitive Icons, colours and grooved triggers ensure suction and irrigation functions are clear.

Minimal Leaks Simplified internal mechanism minimises frustrating leaks.

Configurable Dual shaft attachment options allow for both pistol and trumpet grip.



Single Use Instruments - Scissors

Scissors	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	Metzenbaum Curved Scissors - Standard Blade (Box of 12)	106Y.101	106Y.202	15mm Blade	Double Action Mild Curve Metzenbaum Blade STERILE	1 Box
	Metzenbaum Curved Scissors - Short Blade (Box of 12)	106Y.103	—	10mm Blade	Double Action Curved Metzenbaum Blade STERILE	1 Box










Single Use Instruments - Scissors Insert & Sheath

Scissors	DESCRIPTION	RE-ORDER 33cm	SIZE	FEATURES	MIN QTY
	Metzenbaum Curved Scissors with Insulated Sheath and Locking Nut (Box of 8)	300-600-177S	15mm Blade	Double Action Medium Curved Blade Sheath and Locking Nut included STERILE	1 Box
					




Single Use Instruments - Forceps

Forceps	DESCRIPTION	RE-ORDER 33cm	SIZE	FEATURES	MIN QTY
	Maryland Dissector Forceps with Standard Handle (Box of 12)	107Y.101	19mm Jaw	Double Action Medium Curve Maryland profile STERILE	1 Box
	Maxi Fenestrated Grasper Forceps with Ratchet Handle (Box of 12)	108Y.101	22mm Jaw	Double Action Deep Horizontal Serrations Long Fenestration Rounded Profile STERILE	1 Box
	Johans Fenestrated Grasper Forceps with Ratchet Handle (Box of 12)	108Y.102	17mm Jaw	Double Action Medium Horizontal Serrations Long Fenestration Rounded Profile STERILE	1 Box
	Allis Grasper Forceps with Ratchet Handle (Box of 12)	108Y.103	17mm Jaw	Double Action Small Platform with Deep Vertical Serrations Smooth Inner Jaw Long Tapered Rounded Profile STERILE	1 Box
	Babcock Fenestrated Forceps with Ratchet Handle (Box of 12)	108Y.104	21mm Jaw	Double Action Rounded Tip Platform with Deep Horizontal Serrations Smooth Inner Jaw Tear Drop Fenestration STERILE	1 Box
	Bullnose Forceps with Ratchet Handle (Box of 12)	108Y.105	19mm Jaw	Double Action Deep Horizontal Serrations Traditional Bullnose Rounded Profile STERILE	1 Box
	Upcurve Fenestrated Grasper Forceps with Ratchet Handle (Box of 12)	108Y.108	22mm Jaw	Single Action Deep Horizontal Serrations Double Fenestration Rounded End with Slender Profile STERILE	1 Box





Single Use Instruments - Monopolar Diathermy

Monopolar Diathermy	DESCRIPTION	RE-ORDER 33cm	SIZE	FEATURES	MIN QTY
	L-Hook Electrode (Box of 10)	110Y.101	33cm WL	Standard L-Hook Design Tapered Insulated Tip <div>STERILE</div>	1 Box





Single Use Instruments - Suction Irrigation Devices

Suction Irrigation Devices	DESCRIPTION	RE-ORDER 33cm	RE-ORDER 45cm	SIZE	FEATURES	MIN QTY
	Cadence™ Suction Irrigator - Single Spike (Ctn of 10)	18-003L	18-004L	Ø5mm	Single Spike 3.7m Tubing <div>STERILE</div>	1 Ctn
	Cadence™ Suction Irrigator - Dual Spike (Ctn of 10)	18-001L	18-002L	Ø5mm	Dual Spike 3.7m Tubing <div>STERILE</div>	1 Ctn



Single Use Instruments - Specimen Retrieval Device

Specimen Retrieval Device	DESCRIPTION	RE-ORDER 33cm	SIZE	FEATURES	MIN QTY
	Cadence™ Specimen Retrieval Device (Ctn of 10)	18-005L	200mL Bag	10mm Diameter 29cm Shaft <div>STERILE</div>	1 Ctn
	Cadence™ Specimen Retrieval Device (Ctn of 10)	18-006L	400mL Bag	10mm Diameter 29cm Shaft <div>STERILE</div>	1 Ctn
	Cadence™ Specimen Retrieval Device (Ctn of 10)	18-007L	600mL Bag	10mm Diameter 29cm Shaft <div>STERILE</div>	1 Ctn
	Cadence™ Specimen Retrieval Device (Ctn of 10)	18-008L	1200mL Bag	15mm Diameter 33cm Shaft <div>STERILE</div>	1 Ctn



SUPPORT

With Eisner™ you don't get a call centre, just a dedicated team of specialists ready to tailor your Eisner™ experience.

Experience the difference

Book an obligation free trial

Our team can organise an obligation free, trial assessment of your instruments so you can see and feel the difference before you buy.

Trade-in / replacement program

Manage your tray upgrades

If you are considering decommissioning an instrument our trade-in/replacement service will assist you to address one off replacements or the gradual upgrade of instrument sets. We accept Eisner™ and non-Eisner™ instruments.

Care

Optimise Your Tray Performance

When you choose Eisner™ products our team are available to provide free comprehensive in-service training, support and tray assessments to guide servicing and repair.





Enhance your every move

Warranty | Eisner™ offers a 12 month warranty against manufacturing defect on all metal components of the equipment. Damage and / or breakages caused by improper use or abuse of the equipment are not covered by warranty. Normal wear and tear is not covered. Warranty is null and void where components are mixed with non-Eisner™ instruments.

Service & Repairs | Instruments returned to us for repair/replacement will be assessed and a written quotation provided for you to make a decision on further action. In most cases for instruments beyond repair we can offer a trade-in replacement. Eisner™ keeps factory original spare parts in stock. All repairs should be returned to Eisner™ for repair in order to preserve the integrity of your instruments.

Recommended Lifespan of Instruments | Whilst Eisner™ laparoscopic instruments will function for many years, we recommend replacement prior to any obvious signs of wear and stress which could lead to an issue during surgery. We recommend a lifespan of 3 years on our handles and inserts and suggest that after that time, they are traded in on a new component. Insulated Tubes: By virtue of the fact that insulated tubes will be replaced once the insulation is damaged or compromised in any way, it is not envisioned that the lifespan of a sheath will be beyond 12 months. Instruments should be taken out of circulation immediately if it is found that identifying labelling is no longer visible.



Pacifica Business Park
Unit 18, 97 Canterbury Road
Kilsyth VIC 3137
Australia
T +61 (0)3 9728 4770
F +61 (0)3 9728 6770
eisner@multigate.com.au
www.eisner.com.au



For an obligation
free trial or for more
information call
us today on:
(03) 9728 4770
eisner.com.au

