

# MULTIGATE V-TEX®

## SURGICAL GOWNS

### Uncompromising Safety and Performance.

V-TEX® surgical gowns provide uncompromising performance and the reassurance of compliance with AAMI level 4 standards, making them perfect for long duration, high fluid level procedures.

#### TECHNOLOGY

### Breathable viral barrier.

V-TEX® gowns are made with a unique lightweight tri-laminate material that delivers a high performance impermeable but breathable viral barrier that allows the movement of warm air vapour away from the body whilst stopping bacteria and virus penetration.

#### FEATURES

### AAMI level 4 protection.

Heat sealed seams and a responsive middle layer ensure our V-TEX® gowns meet and exceed AAMI level 4 standards making them the ideal solution for high fluid, long duration procedures.

### Comfort.

With generous dimensions to support upper body movement and a lightweight back panel to assist in regulating the wearer's body temperature, V-TEX® gowns are designed to maximise comfort during procedures.

- ✓ Impervious
- ✓ AAMI level 4 protection
- ✓ Breathable
- ✓ Lightweight back panel
- ✓ Viral barrier
- ✓ Low linting



Lightweight  
and breathable



Protection  
against fluid  
striekthrough  
and viral  
penetration

Standard Gown with 1 Huck Towel				
RE-ORDER	SIZE	GOWN LENGTH	STERILE / NON STERILE	CARTON QUANTITY
27-301	Medium	112cm	Sterile	20 Packs
27-305	Large	125cm	Sterile	20 Packs
27-306	X-Large	135cm	Sterile	20 Packs
27-307	XX-Large	135cm	Sterile	20 Packs
27-302	Large, X-Long	144cm	Sterile	20 Packs
27-303	X-Large, X-Long	148.5cm	Sterile	20 Packs
27-304	XX-Large, X-Long	148.5cm	Sterile	20 Packs

For more information please contact your Business Development Manager:  
**(02) 8718 2888** or visit [multigate.com.au](https://multigate.com.au)