

Our softest, lowest-linting all purpose towel.

A versatile towel that delivers lasting performance especially when wet.

### Features and benefits

Made from non-woven spun lace cellulose/ polyester fabric, MediClean™ is ideal for a wide range of dry or wet CSSD applications. Its high cellulose content makes MediClean™ a favourite with nurses seeking a soft towel that's gentle on their hands. MediClean™ towels are ideal for use as a tray liner and applications where the towel is soaked in water for long periods of time.

- ✓ Lasting performance, even in warm water
- ✓ Exceptional low linting performance
- ✓ Soft and gentle on hands
- ✓ Superior wet and dry strength
- ✓ Absorbs up to 5.5 times its weight in fluid
- ✓ Latex free

#### MediClean™ - Product Specifications (4plys folded)

MATERIAL/ CONSTRUCTION	WEIGHT OF TOWEL	WATER ABSORBENCY	SINKING TIME	TENSILE STRENGTH: Machine Direction	TENSILE STRENGTH: Cross Direction
55% Cellulose 45% Polyester	51gsm	5.5 g/g	<10 sec	97 N	42 N



#### MediClean™ - All Purpose Towel

RE-ORDER	DESCRIPTION	SIZE (CM)	STERILE / NON-STERILE	PIECES PER DISPENSER	PIECES PER CARTON
57-112	MediClean™ Towel Small	33 x 31.8	Non-Sterile	100pcs	1200pcs
57-113	MediClean™ Towel Large	35 x 60	Non-Sterile	100pcs	1200pcs
47-679D	MediClean™ Towel Small (2pcs/pk)	33 x 31.8	Sterile	130pcs (65pks)	520pcs (260pks)
47-680	MediClean™ Towel Large	33 x 31.8	Sterile	-	400pcs
47-681	MediClean™ Towel Large	35 x 60	Sterile	-	250pcs