

# MULTIGATE

## Multi-Strip®



## Wound Closures

### Skin Closure Range

Our Multi-Strip® skin closure adhesive strips are used to keep wound edges together to reduce risk of scarring. It's ideal for securing, closing and supporting surgical incisions, small cuts and lacerations. It assists in increasing the tensile strength of the wound site while providing secure wound support. It offers a better cosmetic outcome than sutures or staples for superficial wounds.



- ✓ Sterile & breathable microporous material
- ✓ Non-woven, reinforced strips for added strength
- ✓ Easy application and removal
- ✓ Hypoallergenic adhesive
- ✓ Can be used together with skin staples or post removal to support wound healing
- ✓ Allows finer wound edge approximation
- ✓ Latex free

#### Multi-Strip® Skin Closures

| RE-ORDER | DESCRIPTION                            | SIZE         | STERILE / NON-STERILE | PACK PER BOX | BOX PER CARTON | PACK PER CARTON |
|----------|--|--------------|-----------------------|--------------|----------------|-----------------|
| 30-351   | Multi-Strip® Skin Closure (Pack of 5 ) | 3mm x 75mm   | Sterile               | 50           | 4              | 200pks          |
| 30-353   | Multi-Strip® Skin Closure (Pack of 6)  | 6mm x 38cm   | Sterile               | 50           | 4              | 200pks          |
| 30-352   | Multi-Strip® Skin Closure (Pack of 3)  | 6mm x 7 5mm  | Sterile               | 50           | 4              | 200pks          |
| 30-354   | Multi-Strip® Skin Closure (Pack of 10) | 6mm x 100mm  | Sterile               | 50           | 4              | 200pks          |
| 30-355   | Multi-Strip® Skin Closure (Pack of 6)  | 13mm x 100mm | Sterile               | 50           | 4              | 200pks          |
| 30-356   | Multi-Strip® Skin Closure (Pack of 4)  | 25mm x 100mm | Sterile               | 25           | 4              | 100pks          |