

# Needle Counters

## Improving clinician safety.

When we saw the challenges you face with sharps management we created a range of needle counters with your safety in mind. Our new slide lock needle counters are designed to reduce the risk and incidence of injury during sharps management.

### DESIGN

## Safe Handling and Disposal of Sharps.

A simple slide lock reduces the risk of injury from accidental opening during the disposal process and the option of a clear lid enables an effective sharps count with the lid closed.

- ✓ Safe disposal
- ✓ Reduced injury risk
- ✓ Firm instrument grip
- ✓ Safe visual count
- ✓ Ergonomic design
- ✓ Recyclable plastic



### FEATURES

## High Density Foam.

The high density foam construction of our counters resists distortion from the sterilisation process and provides a stronger grip on sharps, allowing for a safer and more efficient work flow.

## Tailored to You.

Our needle counters are available in various combinations of magnetic and foam styles and come in a range of different counts.


Our counters can be included in Custom Procedure Packs or purchased as a sterile single, ensuring we have a solution to suit your specific needs.

## Recyclable Plastic.

The entire range of needle counters are made from recyclable plastic.




Magnetic	
RE-ORDER Sterile	<b>19-231</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-231NS</b>
CARTON QTY	150 Pieces




Magnetic 30 Count Foam Strip Clear Lid	
RE-ORDER Sterile	<b>19-237</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-237NS</b>
CARTON QTY	200 Pieces




Magnetic 15 Count Foam Strips	
RE-ORDER Sterile	<b>19-232</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-232NS</b>
CARTON QTY	200 Pieces




Magnetic 30 Count Foam Block	
RE-ORDER Sterile	<b>19-238</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-238NS</b>
CARTON QTY	200 Pieces




Magnetic 20 Count Foam Strip Clear Lid	
RE-ORDER Sterile	<b>19-233</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-233NS</b>
CARTON QTY	200 Pieces




Magnetic 40 Count Foam Block	
RE-ORDER Sterile	<b>19-239</b>
PER BOX	10 Packs
PER CARTON	4 Boxes
CARTON QTY	40 Packs
RE-ORDER Non-sterile	<b>19-239NS</b>
CARTON QTY	100 Pieces



Magnetic 20 Count Foam Block	
RE-ORDER Sterile	<b>19-234</b>
PER BOX	10 Packs
PER CARTON	8 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-234NS</b>
CARTON QTY	100 Pieces



Magnetic 42 Count Foam Strip	
RE-ORDER Sterile	<b>19-240</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-240NS</b>
CARTON QTY	200 Pieces



Magnetic 30 Count	
RE-ORDER Sterile	<b>19-236</b>
PER BOX	20 Packs
PER CARTON	4 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-236NS</b>
CARTON QTY	200 Pieces



Magnetic 20 Count Foam Strip Clear Lid	
RE-ORDER Sterile	<b>19-241</b>
PER BOX	10 Packs
PER CARTON	8 Boxes
CARTON QTY	80 Packs
RE-ORDER Non-sterile	<b>19-241NS</b>
CARTON QTY	100 Pieces