



Protective Eyewear

DESCRIPTION

Protective eyewear is a non-sterile device intended to be worn over the eyes for protection from blood and other body fluid

splashes or foreign matter during surgical, clinical or laboratory procedures.

54-100NS (BLACK FRAME)



54-106NS (CLEAR FRAME)



FEATURES AND BENEFITS

- ✓ AS1337.1-2010 compliant
- ✓ High quality black & extendable frame
- ✓ Optically transparent lenses
- ✓ Anti-scratch lens
- ✓ Medium impact protection, suitable for outdoor use
- ✓ Light weight & Ergonomically designed
- ✓ Wall mounted dispenser box

- ✓ AS1337.1-2010 compliant
- ✓ Shaped side arms for fitting over prescription eyewear
- ✓ One-piece lens
- ✓ Clear anti-fog & anti-scratch lens
- ✓ Suitable for outdoor use
- ✓ Vented temple

PRODUCT DATA

- ARTG ID 235988
- Transparent polycarbonate lens
- Matt black nylon frame & temple
- Nose pad: Thermoplastic Rubber (TPR)
- Latex free

- ARTG ID 235988
- Transparent polycarbonate lens
- Transparent polycarbonate side shield
- Transparent polycarbonate nose pad

Personal Protective Equipments

RE-ORDER	DESCRIPTION	STERILE / NON-STERILE	PIECES PER INNER BOX	INNER PER CARTON	TOTAL PIECES PER CARTON
54-100NS	Protective Eyewear (Black Frame)	Non-Sterile	15pcs	10 Boxes	150pcs
54-106NS	Protective Eyewear (Clear Frame)	Non-Sterile	12pcs	25 Boxes	300pcs