

MULTIGATE SafeCount™



Needle Counter 20 count / 20 capacity

STERILE

PRECISION,
ORGANISATION
AND PROTECTION

The SafeCount™ needle counter securely manages and counts sharps with innovative features designed for precision, organisation, and enhanced protection.

RE-ORDER CODE

CEA-0014S

FEATURES AND BENEFITS

- ✓ Clear lid allows for sharps to be viewed and counted whilst the needle counter is closed
- ✓ Printed numbers to assist with counting sharps
- ✓ X-ray detectable secure lock ensures needle counter is closed and ready for safe disposal
- ✓ Adhesive fixation on rear to keep needle counter in place
- ✓ Easy to identify product is yellow, the industry standard colour for sharps



PRODUCT DETAILS

Description	Needle counter 20 count / 20 capacity Hybrid with clear lid in shallow option	Size when closed	16 x 127 x 74mm
Material	Injection moulded polystyrene, magnet sheet, PE foam and paper	Units per box	20
Latex free	Yes	Boxes per carton	4
Shelf life	—	Carton quantity	80
Non-Sterile code	CEA-0014	Sterile	Sterile
		Use	Single use
		Colour	Yellow, universal colour to indicate sharps



As part of our commitment to better product stewardship, all SafeCount™ needle counters sold to Australian customers are supplied Climate Neutral¹.

Disclaimers: 1. Multigate only purchases carbon credits via Carbon Neutral in approved, registered projects which have been independently assessed and verified as creating additional resources that help combat climate change using avoidance or sequestration means.

IP Protected Product: Intellectual Property rights such as patents, patent applications, registered designs, trademarks, etc. apply to this product group. See website for details.

FOR MORE INFORMATION

multigate.com.au or call 1800 023 420 (Australia)

multigate.co.nz or call 0800 880 013 (New Zealand)