

MULTIGATE

SoftPro®

SURGICAL GOWNS



Advanced Protection.

SoftPro® surgical gowns offer superior protection and comfort. Their 5-layer polypropylene construction delivers an effective splash resistant barrier against bacterial and fluid transfer whilst maintaining breathability, ensuring you stay comfortable and cooler during OR procedures.

TECHNOLOGY

Protection.

SoftPro® gowns are manufactured with a high quality five-layer SSMMS fabric that delivers a high-performance splash resistant barrier that repels fluids and bacteria.

FEATURES

Comfort.

The SoftPro® fabric is low linting and easily draped. With generous sleeve sizes and range of lengths, the SoftPro® range has been designed with your comfort in mind.

AAMI level 3 protection.

SoftPro® gowns meet and exceed AAMI level 3 standards making them the ideal solution for procedures involving light to moderate fluid levels.

Reinforced.

SoftPro® reinforced gowns feature a strong, yet breathable extra layer of reinforcement in the critical arm and front areas for added protection against fluid strikethrough.

- ✓ Comfortable
- ✓ Breathable
- ✓ Splash resistant barrier
- ✓ AAMI level 3 protection
- ✓ Reinforced critical areas
- ✓ Low linting



Microbial barrier that delivers improved barrier performance



Lightweight & breathable

Standard Gowns

RE-ORDER Gown +2 Towels	RE-ORDER Gown +1 Huck Towel	SIZE	STERILE / NON-STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
98-001	98-001H	Medium	Sterile	20 Packs	99-101
98-002	98-002H	Large	Sterile	20 Packs	99-102
98-003	98-003H	X-Large	Sterile	20 Packs	99-103
98-004	98-004H	XX-Large	Sterile	20 Packs	99-104
98-006	98-006H	XXX-Large	Sterile	20 Packs	99-106
98-005	-	X-Large, X-Long	Sterile	20 Packs	99-105
98-022	66-413H	Large *2 gown pack	Sterile	10 Packs	99-102

*Includes two large gowns and either 4 standard towels or 2 Huck Towels

Reinforced Gowns

RE-ORDER Gown +2 Towels	RE-ORDER Gown +1 Huck Towel	SIZE	STERILE / NON-STERILE	PACKS PER CARTON	RE-ORDER Non-sterile Gown only
99-001	-	Medium	Sterile	20 Packs	-
99-002	99-002H	Large	Sterile	20 Packs	99-202
99-003	99-003H	X-Large	Sterile	20 Packs	99-203
99-004	99-004H	XX-Large	Sterile	20 Packs	99-204
99-006	-	XXX-Large	Sterile	20 Packs	99-206

For more information please contact your Business Development Manager:

(02) 8718 2888 or visit multigate.com.au