

Sterilisation Apparel & Accessories



CSSD Accessories and Consumables

At Multigate we combine our deep understanding of how you work with our design capabilities to innovate and better serve you, your patients and our partners.

For 30 years we've created medical products, services and experiences that give you the advantage to perform at your best. We continue to do this.

Inspiring performances and better patient outcomes that will lead the world of hospital care to new standards.

We are Multigate.

We understand that sterile management, clinician protection and streamlined daily workflows are paramount in the CSSD.

With that in mind, we carefully designed our range of protective apparel, towels and tray liners and sterilisation accessories to power your performance.



Protective Apparel and Eyewear

- Thumb Secure Isolation/Protective Gown
- Liteback Thumb Secure Protective Gown
- Protective Eyewear

Thumb Secure Isolation/Protective Gown Non-Sterile					
RE-ORDER	COLOUR	SIZE	PIECES/BOXES	BOXES/CARTON	TOTAL CARTON QTY (PIECES)
29-285	Yellow	Regular	15	5	75
29-286	Yellow	X-Large	15	5	75
29-385	Blue	Regular	15	5	75
29-386	Blue	X-Large	15	5	75

Liteback Thumb Secure Protective Gown Non-Sterile					
RE-ORDER	COLOUR	SIZE	PIECES/BOXES	BOXES/CARTON	TOTAL CARTON QTY (PIECES)
29-291	Yellow	Regular	10	6	60
29-293	Yellow	XX-Large	10	6	60

Protective Eyewear Non-Sterile				
RE-ORDER	COLOUR	PIECES/BOXES	BOXES/CARTON	TOTAL CARTON QTY (PAIRS)
54-100NS	Black	15	10	150



Thumb Secure Protective Gown



Liteback Thumb Secure Protective Gown



Protective Eyewear

Towels and Tray Liners

- MediClean
- MediCel®
- All Purpose Towels

MediClean Non-Sterile			
RE-ORDER	SIZE	PIECES/BOXES	BOXES/CARTON
57-112	Small	100	12
57-113	Large	100	12

MediCel® Non-Sterile			
RE-ORDER	SIZE	PIECES/BOXES	BOXES/CARTON
57-120	Small	100	12
57-121	Large	100	12

All Purpose Towels Non-Sterile			
RE-ORDER	SIZE	PIECES/BOXES	BOXES/CARTON
57-885	Small	100	12
57-886	Large	100	12



MediClean



MediCel®




All Purpose Towels

Pro-Zone® Sterilisation Wraps

Versatile high quality wrapping system that delivers many performance advantages including:

- Strength and resistance to punctures, tears and abrasions*
- Anti-static coating reduces static “shocks” to clinicians
- Low linting
- Soft and easy to fold
- Excellent drapability and low memory
- Excellent microbial barrier
- Effective sterilant penetration

* when correct handling techniques are used



Pro-Zone®
The original Pro-Zone® sterilisation wrap.

- Weights:**
 - 47gsm
 - 61gsm
 - 75 gsm
- Compatible with:**
 - Steam
 - ETO
 - Low temperature gas plasma
 - Low temperature steam/formaldehyde
- Standards Compliance:**
 - AS/NZS 4187:2014
 - ISO 11607



Pro-Zone®+2
Pro-Zone®+2 boasts the innovative use of 2 coloured sheets – blue and lilac, bonded together, enabling clinicians to rapidly identify any tears or damage in the wrap.

- Weights:**
 - 47gsm
 - 61gsm
 - 75 gsm
- Compatible with:**
 - Steam
 - ETO
 - Low temperature gas plasma
 - Low temperature steam/formaldehyde
- Standards Compliance:**
 - AS/NZS 4187:2014
 - ISO 11607



Pro-Zone®+4G
Pro-Zone®+4G offers the protection of 2 bonded sheets, plus the added advantage of an absorbent layer to wick any residual condensate that may compromise the integrity of the wrap.


- Weights:**
 - 61gsm
 - 75 gsm
- Compatible with:**
 - Steam
 - ETO
 - Low temperature steam/formaldehyde
- Standards Compliance:**
 - AS/NZS 4187:2014
 - ISO 11607

Pro-Zone® Accessories

- Protective Corner
- Transfer Tray

Protective Corner Non-Sterile		
RE-ORDER	COLOUR	CARTON QTY (PIECES)
DP 2704	Orange	300

Transfer Tray Non-Sterile			
RE-ORDER	SIZE (CM)	COLOUR	CARTON QTY (PIECES)
DP 2702	40 x 30.3 (small)	Green	10
DP 1702	56 x 30.3 (large)	Green	10



Protective Corner



Transfer Tray

Empowering
performance
since 1986

Multigate Medical Products

27 Llewellyn Avenue

Villawood NSW 2163

Australia

T +61 (0)2 8718 2888

F +61 (0)2 8718 2889

enquiries@multigate.com.au

www.multigate.com.au

For more information please contact your
local Business Development Manager:

(02) 8718 2888 or visit **multigate.com.au**