



Empowering performance since 1986

Multigate Medical Products
 27 Llewellyn Avenue
 Villawood NSW 2163
 Australia
 T +61 (0)2 8718 2888
 F +61 (0)2 8718 2889
 enquiries@multigate.com.au
 www.multigate.com.au

Trauma/Burns Drape



Featuring a 270 Degree Fluid Control Pouch
Innovation Focus: Fluid Management

v2.0

Your expert opinion matters
 Your feedback will help us to improve our products and power your performance. We would appreciate a quick appraisal.

What do you think of this product on a scale of 1-5 (1 being very poor and 5 being highly effective)

Easy to use?	Performance enhancing?	Would you recommend?
<input type="radio"/> 1	<input type="radio"/> 1	<input type="radio"/> Yes
<input type="radio"/> 2	<input type="radio"/> 2	<input type="radio"/> Maybe
<input type="radio"/> 3	<input type="radio"/> 3	<input type="radio"/> No
<input type="radio"/> 4	<input type="radio"/> 4	
<input type="radio"/> 5	<input type="radio"/> 5	

For more information please contact your local Business Development Manager:
(02) 8718 2888



The above image is a suggestion on how the drapes can be combined for a patient who requires full body treatment. However, each drape can be used separately depending on what is required for the individual case presented.

A drape that puts you in control of large volumes of fluid in trauma and burns cases.



01 | EFFECTIVE FLUID MANAGEMENT

Our unique 270-degree fluid control pouch design enables you to effectively contain large volumes of fluid to maintain a safe & sterile environment and to assist in monitoring the amount of fluid & blood loss.



02 | IMPERVIOUS

The entire drape is impervious, preventing fluid strikethrough and assisting in the management of any fluid overflow by channelling the fluid into the surrounding pouch.



03 | AAMI LEVEL 4 PROTECTION

Constructed using our state of the art Protex@4 fabric, our drapes provide AAMI level 4 protection across the entire drape (not just in the critical zones), delivering maximum clinician and patient protection against strikethrough and microorganism penetration.



04 | SUPERIOR ABSORBENCY

Our trauma/burns drapes support effective fluid management by being absorbent across the entire drape, with super absorbency in the critical zones.

- ✓ 270-degree pouch design to manage large volumes of fluid
- ✓ Entire drape is impervious
- ✓ AAMI level 4 protection
- ✓ Entire drape is absorbent with extra absorbency in the critical zones

Trauma/Burns Drape

RE-ORDER	DESCRIPTION	SIZE	CARTON QUANTITY	NON-STERILE
35-168	Trauma/Burns Drape with U-Split	230cm x 270cm	12	35-168NS
35-177	Trauma/Burns Drape - Small	200cm x 200cm	20	35-177NS



Application mitts enable the accurate application of the drape whilst still maintaining sterility.

Multiple suction ports to assist with efficient drainage of the pouch when required.