

Trauma/Burns Drape

A drape that puts you in control of large volumes of fluid in trauma and burns cases.

DESIGN

EFFECTIVE FLUID MANAGEMENT

Our unique 270-degree fluid control pouch design enables you to effectively contain large volumes of fluid to maintain a safe & sterile environment and to assist in monitoring the amount of fluid & blood loss.

IMPERVIOUS

The entire drape is impervious, preventing fluid strikethrough and assisting in the management of any fluid overflow by channelling the fluid into the surrounding pouch.

AAMI LEVEL 4 PROTECTION

Constructed using our state of the art Protex®4 fabric, our drapes provide AAMI level 4 protection across the entire drape (not just in the critical zones), delivering maximum clinician and patient protection against strikethrough and microorganism penetration.

SUPERIOR ABSORBENCY

Our trauma/burns drapes support effective fluid management by being absorbent across the entire drape, with super absorbency in the critical zones.

FEATURES

- ✓ 270-degree pouch design to manage large volumes of fluid
- ✓ Entire drape is impervious
- ✓ AAMI level 4 protection
- ✓ Entire drape is absorbent with extra absorbency in the critical zones



Application mitts enable the accurate application of the drape whilst still maintaining sterility.



Multiple suction ports to assist with efficient drainage of the pouch when required.



STERILE RE-ORDER	DESCRIPTION	SIZE	CARTON QUANTITY	NON-STERILE RE-ORDER
35-168	Trauma/Burns Drape with U-Split	230cm x 270cm	12 Packs	35-168NS
35-177	Trauma/Burns Drape - Small	200cm x 200cm	20 Packs	35-177NS
35-126	Basic Trauma Drape Sheet	150cm x 80cm	40Packs	35-126NS