

Under-pads

Absorbing the pain of unwanted fluids.

Addressing nurses concerns over fluid management we created an Under-pad range designed to maximise absorbency and reduce unwanted leaks.



DESIGN

The Multigate range of under-pads are designed to maximise fluid wicking and prevent leaks.

FEATURES & BENEFITS

- ✓ The top non-woven layer is soft and non-irritating for maximum comfort.
- ✓ Five-ply and six-ply design available to meet your different absorbency management requirements.
- ✓ 4 sealed sides prevent unwanted leaks and contamination.
- ✓ Unique design of embossed honeycomb pattern provides faster fluid absorption.
- ✓ A waterproof backing to prevent strikethrough
- ✓ 8% more fluid holding Capacity *

Unique
Honeycomb
Pattern

Multigate Under-Pad Non-sterile

RE-ORDER	DESCRIPTION	SIZE (CM)	PIECES/ POLYBAG	PIECES/ CARTON
22-121	Under-pad: 5 ply	55 x 40	50	300
22-124	Absorbent Pad	40 x 30	50	600
22-128	Under-pad Plus: 6 ply	55 x 40	50	300
22-129	Under-pad Plus: 6 ply	40 x 30	50	600

* Internal testing data, compared with other traditional equivalent bluey products.