

Under-pads

Absorbing the pain of unwanted fluids.


Addressing nurses concerns over fluid management we created an Under-pad range designed to maximise absorbency and reduce unwanted leaks.

DESIGN

With their unique five and six-layer construction, Multigate Under-pads are designed to maximise fluid wicking, prevent leaks and minimise the impact of unwanted fluids.

Under-Pad Non-sterile

RE-ORDER Sterile	DESCRIPTION	SIZE (CM)	PIECES/ POLYBAG	PIECES/ CARTON
22-121	Under-pad: 5 ply	55 x 40	50	300
22-124	Absorbent Pad	40 x 30	50	600
22-128	Under-pad Plus: 6 ply	55 x 40	50	300



FEATURES

Highly absorbent.

Our five-ply and six-ply Under-pads are made from a soft, non-woven fabric and a highly absorbent fill material to maximise their ability to contain unwanted fluids.

Leak prevention.

Sealed on four sides, with a waterproof backing to prevent strikethrough, Multigate Under-pads are designed to avoid the added pain of leaks.

