

Dressing Pack

08-559NP (Previous Version)



Basic Dressing Pack with 5 Non-Woven Swabs (Tear Pack)

Re-order code: **08-559NP** (Previous Version)

| | |
|-------------------|------------|
| Size | - |
| Units per Inner | 20pks |
| Inners per Carton | 8 bags |
| Carton Quantity | 160pks |
| Sterile | Sterile |
| Use | Single use |
| Colour | - |

Pack Contents

| Description | Size | Quantity in Pack |
|------------------------|--------------------------|------------------|
| Sterile Field | 49cm x 49cm | 1 |
| Dressing Tray Blue | 13.7cm x 13.7cm x 1.63cm | 1 |
| Paper Towel | 40cm x 38cm | 1 |
| Non-Woven Swabs | 7.5cm x 7.5cm | 5 |
| Plastic Forceps Yellow | 10.8cm | 1 |
| Plastic Forceps Blue | 10.8cm | 2 |