


Needle Counter

19-241 (Previous Version)



Needle Counter - 20 Count,
Foam Strip, Clear Lid

Re-order code: 19-241 (Previous Version)

Size	-
Units per Inner	10 pks
Inners per Carton	8 boxes
Carton Quantity	80 pks
Sterile	Sterile
Use	Single use
Colour	-

DESIGN

Safe Handling and Disposal of Sharps.

A simple slide lock reduces the risk of injury from accidental opening during the disposal process and the option of a clear lid enables an effective sharps count with the lid closed.

For added patient safety, the slide lock is also x-ray detectable for post-operative detection.

- ✓ Safe disposal
- ✓ Reduced injury risk
- ✓ Firm instrument grip
- ✓ Safe visual count
- ✓ Ergonomic design
- ✓ Recyclable plastic