

# MULTIGATE

## Protex<sup>®</sup>4



Advanced Protection System  
**Innovation Focus: Barrier Performance**

V2.0



At Multigate we combine our deep understanding of how you work with our design capabilities to innovate and better serve you, your patients and our partners.

For over 30 years we've created medical products, services and experiences that give you the advantage to perform at your best. We continue to do this.

Inspiring performances and better patient outcomes that will lead the world of hospital care to new standards.

We are Multigate.



## ADVANCED PROTECTION SYSTEM

We have combined over fifteen years of experience in drape development for the Australian perioperative market, with state of the art materials and manufacturing technologies, to bring you the latest in imperviousness and absorbency, delivering unbeatable patient and surgical team protection.


- ✓ Impervious
- ✓ High absorption capacity and speed
- ✓ Soft and pliable
- ✓ Strong and durable
- ✓ Low linting
- ✓ No latex



Improved microbial barrier performance



Excellent draping capability




## TECHNOLOGY — STATE OF THE ART MATERIALS TECHNOLOGY

Multigate Protex<sup>®</sup>4 drapes are made from a pliable yet tough impervious bi-laminate material supported by an additional ultra-reinforcement layer to further boost absorbency if required. They are designed to eliminate bacterial transfer and deliver maximum absorbency and strikethrough prevention.

## FEATURES — SUPERIOR BARRIER PERFORMANCE

Superior Barrier Performance Multigate Protex<sup>®</sup>4 drapes deliver impervious AAMI level 4 protection not just in the critical zone but across the entire drape.



## High performance tethers.

Our tethers provide a strong heavy duty hold to manage quivers and leads with confidence.

## Reliable adhesion.

Protex<sup>®</sup>4 drapes offer a reliable incise film and adhesion performance and are suitable for all skin types.

## Fluid collection.

Fluid pouches made from polyethylene offer a durable heavy duty seal and are combined with strategically placed suction ports to aid the capture, monitoring and management of fluid and blood loss.

## For more information on Protex<sup>®</sup>4 and our drapes range



visit [multigate.com.au/protex-4](https://www.multigate.com.au/protex-4)

Alternatively, contact your local Business Development Manager on (02) 8718 2888.

[www.multigate.com/qr05](https://www.multigate.com/qr05)

Empowering  
performance  
since 1986

**Multigate Medical Products**

27 Llewellyn Avenue

Villawood NSW 2163

Australia

T +61 (0)2 8718 2888

F +61 (0)2 8718 2889

[enquiries@multigate.com.au](mailto:enquiries@multigate.com.au)

[www.multigate.com.au](http://www.multigate.com.au)

For more information please contact your  
local Business Development Manager:

**(02) 8718 2888** or visit [multigate.com.au](http://multigate.com.au)